Harvard MTA-PT OptiCaps®

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Harvard MTA-PT OptiCaps® is an endodontic repair cement in OptiCaps®. Due to its unique flowable consistency in combination with a short setting time Harvard MTA-PT is the ideal choice for **pulpotomy**.

Harvard MTA-PT powder consists of very fine hydrophilic particles of several hydraulic mineral oxides. After contact with the liquid it forms a gel that hardens to an impermeable barrier.

Harvard MTA-PT OptiCaps® are easily activated and the content of the capsule is easily ejected with the Harvard Applier OptiCaps®. Mixing of Harvard MTA-PT OptiCaps® (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/minute.

Primary indication:

- Pulpotomy

Secondary indication:

- Pulp capping (direct and indirect)

Contraindications / side-effects

None known.

1. Activation and mixing

See "Instructions for the activation and mixing of Harvard OptiCaps®"

Mixing time for the Harvard MTA-PT OptiCaps® is 30 seconds.

Attention:

Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting and lag times may impair or prevent application of the material.

The material must be extruded within **10 seconds** after the end of mixing.

To prevent dehydration during setting, apply Harvard MTA-PT immediately after mixing.

2. Working time

Working time at 23° C (from the start of mixing):
about 2 minutes

3. Applications**3.1. Pulpotomy**

Place rubber dam. After complete removal of all carious infected areas remove the infected coronal pulp. Stop any bleeding that may occur

Mix a capsule of Harvard MTA-PT OptiCaps® as described under item 1 and extrude it on a glass plate.

Using a suitable instrument cover the entire surface of the exposed pulp with Harvard MTA-PT.

At the earliest **3 minutes** after application of Harvard MTA-PT proceed with tooth restoration as usual.

Check pulp vitality and status regularly.

3.2. Pulp capping

Place rubber dam and prepare the cavity. Rinse the cavity and exposed pulpal areas with a suitable disinfectant.

Mix a capsule of Harvard MTA-PT OptiCaps® as described under item 1 and extrude it on a glass plate.

With a suitable instrument apply a small amount of Harvard MTA-PT over the exposed pulp and remove excess moisture with a cotton pellet.

At the earliest **3 minutes** after application of Harvard MTA-PT place a small amount of a flowable light cure liner (e.g. Harvard IonoLine) and light cure.

Etch the remaining cavity walls according to the total-etch-technique with Harvard Etch and apply a suitable bonding agent (e.g. Harvard Bond TE Bond) according to the corresponding instructions.

Place a light cure composite (e.g. Harvard UltraFill) according to the instructions and light cure.

Check pulp vitality and status regularly.

Notes

- In the first hour after application handle the placed MTA cement carefully.
- Store Harvard MTA-PT OptiCaps® in the sealed packaging in a dry place prior to use.
- Harvard MTA-PT can cause discoloration.

Storage

Store Harvard MTA-PT at a dry place at 2 - 25 °C (36 °F - 77 °F). Do not use after expiry date.

Harvard MTA-PT OptiCaps® are for single use only.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!**Store product out of reach of children!**

Keep the instructions for use for the duration of the application.

Order-No. Article

7081512 Harvard MTA-PT 2 OptiCaps® ea. 0.25 g, ea. packed in an aluminium pouch

7092000 Harvard Applier OptiCaps®

Instructions for the activation and mixing of Harvard MTA-PT OptiCaps®

	<p>Harvard OptiCaps® before activation.</p>
	<p>For activation of the Harvard OptiCaps®, press the plunger on a hard and plane surface to the end into the capsule.</p>
<p>Click 1x</p>	<p>Insert the Harvard OptiCaps® into the Harvard Applier OptiCaps® and click once. Note: The plunger must be at the same level as the bottom of the capsule.</p>
	<p>Insert the Harvard OptiCaps® into a mixer, close lid and mix immediately for 30 seconds at 4,300 oscillations/minute.</p>
<p>1</p> <p>Click 2x</p> <p>2</p>	<p>Insert the Harvard OptiCaps® into the Harvard Applier OptiCaps®. Remove the pin from the nozzle. If you leave it, the capsule can burst.</p> <p>Pull the lever 2 times (2 clicks) to prime the the Harvard OptiCaps®.</p>
<p>A</p>	<p>Extrude the mixed material on a glass plate or apply directly. Unlock the Harvard Applier OptiCaps® (press button A) and remove the Harvard OptiCaps®.</p> <p>The optimal amount of mixed material is guaranteed only with the Harvard Applier OptiCaps® (Order-No. 7092000).</p>

For the selection of a suitable capsule mixer, our sales and marketing colleagues are gladly available to you.