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7083607 Harvard Bond TE Mono, 5ml Flasche

Instructions for Use EN
Light curing total-etch bonding

Harvard Bond TE Mono is a simple to use light cure one component bonding agent for the total etch technique with a phosphoric acid etching gel ("Etch and Rinse"). It is designed for strong bonding of light cure composites to etched enamel and dentine and to nonprecious and precious metals.

Good, long lasting adhesive strength and good biocompatibility are attained by polycarboxylic acids.

Harvard Bond TE Mono is compatible with all current brands of visible light cure methacrylate-based composite materials.

It is ethanol-based and hydrophilic. Harvard Bond TE Mono can be used on slightly moist dentin surfaces (wet bonding technique).

Indications

Harvard Bond TE Mono is designed for bonding of light cure composites to enamel and dentine

Contraindications

The placement of Harvard Bond TE Mono is contraindicated if the recommended application technique is not possible. Also do not use Harvard Bond TE Mono if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Harvard Bond TE Mono may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the material. Neither store the material in proximity of eugenol containing products, nor let the material come into contact with materials containing eugenol.

Application
1. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. Margins should have a slight (0.5-1.0 mm) bevel placed in the enamel to increase the surface area for greater bond strength.

2. Pulp protection

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard CalciumHydroxide).

3. Application of Harvard Bond TE Mono
3.1. Enamel and Dentine Conditioning

According to the total etch technique apply a 37% phosphoric acid etching gel (e.g. Harvard Etch) onto the enamel and dentine surfaces beginning with the enamel bevels. Condition the enamel for at least 15 seconds and the dentine for 15 seconds or less. Deciduous teeth are etched correspondingly longer. Rinse for 20 seconds with water. Dry it in a water-free and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of Harvard Bond TE Mono. The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

3.2. Application of Harvard Bond TE Mono

Apply Harvard Bond TE Mono generously with a brush onto the enamel and dentin surfaces for 30 seconds with agitation. The material should build a homogeneous layer. Remove excess material carefully. Dry cautiously with oil free air for about 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. Do no desiccate the dentine.

Note:

If not used immediately, place dispensed Harvard Bond TE Mono in subdued light to prevent premature polymerization by incident light. Harvard Bond TE Mono will not selfcure.

Then light cure with a suitable polymerization unit (wavelength range 400-500 nm, light intensity of at least 500mW/cm²) for 20 seconds before application of a second layer of Harvard Bond TE Mono.

3.3. Application of a Second Layer of Harvard Bond TE Mono

Apply again the Harvard Bond TE Mono generously with a brush onto the adhesive surfaces as described under point 3.2. before placement of a light cure composite. **Do not touch the treated surface!**

4. Placement of the restorative

Apply the restorative material according to the instructions of the manufacturer. Best results are obtained with application of a thin layer of a light cure flowable composite (e.g. Harvard RestoreFlow) followed by the application of a moldable composite (e.g. Harvard Restore) according to the instructions.

5. Warnings

- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard Bond TE Mono may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!
Store product out of reach of children!

Keep the instructions for use for the duration of the application.

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