## Harvard InterLock

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## Harvard InterLock

### Light-curing self etching adhesive, Primer & Bond

Harvard InterLock is a simple to use self-etching light cure bonding system. It consists of a self-etching primer Harvard InterLock Prime and the bonding agent Harvard InterLock Bond. The application of the Harvard InterLock components one after the other guarantees safe and reproducible results without technique sensitivity.

#### Indications

Harvard InterLock is designed for bonding of light cure composites to:

- enamel and dentine
- silicate ceramics
- zirconia
- non-precious metals

#### Contraindications

The placement of Harvard InterLock is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard InterLock if the patient is known to be allergic to any of the ingredients.

#### Side effects

In singular cases, Harvard InterLock may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calcium hydroxide material.

#### Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol can impair the hardening of the adhesive. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

### Application

#### 1. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

Margins should have a slight (**0.5 - 1.0 mm**) bevel placed in the enamel to increase the surface area for greater bond strength.

#### 2. Pulp protection

Cavity floor of deep excavations should be covered with a thin layer of calcium hydroxide material.

#### 3. Application of Harvard InterLock

##### 3.1 Harvard InterLock Prime

Apply Harvard InterLock Prime with the supplied applicator onto the **slightly wet** enamel- and dentin surfaces (as they are usual in clinical procedures) for about 10 seconds with agitation.

The material should build a homogeneous layer. Remove excess material carefully. Air thin gently to remove the volatile components and to disperse the adhesive. Take care, that the adhesive is not removed from the tooth surface by a strong air stream.

**Do not touch the treated surface!**

##### 3.2. Harvard InterLock Bond

Apply Harvard InterLock Bond for about 10 seconds with agitation using the applicator. Air thin to remove the volatile components and to disperse the bonding to an even layer until the dentine surface appears slightly dull (the enamel surface remains glossy). Then light cure with a suitable polymerization unit (wavelength range 400–500 nm, light intensity of at least 500mW/cm<sup>2</sup>) for 20 seconds before placement of a light cure composite.

#### Note:

Application of Harvard InterLock Prime and Harvard InterLock Bond can be done with the same applicator.

#### 4. placement of the restorative

Apply the restorative material according to the instructions of the manufacturer.

Best results are obtained with application of a thin layer of a light cure flowable composite (e.g. Harvard PremiumFlow) followed by the application of a moldable composite (e.g. Harvard PremiumFill) according to the instructions.

#### 5. Warnings

- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard InterLock may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

#### Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

#### Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Harvard Dental Internationals sole obligation shall be repair or replacement of the Harvard Dental International product.

#### Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### Only for dental use!

#### Store product out of reach of children!

The instructions for use has to be kept for the duration of the application.

#### Order-No. Article

Harvard InterLock  
7083603 2 x 5 ml bottle, 50 micro applicators,  
5 mixing wells