Ionomer Liner, lighthärtend

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Ionomer liner, light-cure

Harvard IonoLine is a light cure ionomer liner. The good radiopacity of the material ensures good postoperative diagnosis. Harvard IonoLine releases fluoride.

Indications

Lining of cavities for following filling procedures with light cure composites.

Contraindications

None known.

Side effects

In singular cases, Harvard IonoLine may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor allow not polymerized composite to come into contact with materials containing eugenol.

Application

1. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction.

For pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material (e.g. Harvard CalciumHydroxide).

Apply a suitable bonding agent (e.g. Harvard Bond SE Mono) to dentine and enamel.

2. Application of Harvard IonoLine

Replace the cap of the syringe by a Harvard NeedleTip H20.

Note:

Squeeze out the first small drop on a pad to ensure proper handling.

Apply Harvard IonoLine in a thin layer until the desired thickness is reached (**max. 2 mm**). Do not apply to the cavity edges and remove any paste sticking to them.

Light cure Harvard IonoLine with a suitable polymerization unit (wavelength range 400–500 nm, light intensity of at least 500mW/cm²) for **20 seconds**.

Hold the light emission window as close as possible to the filling material.

Proceed with placement of a light cure composite (e.g. Harvard PremiumFill) according to the corresponding instructions.

3. Warnings

- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard IonoLine may have an irritant effect and can lead to sensitization against methacrylates.

- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Close syringe immediately after use to avoid exposure to light.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Keep the instructions for use for the duration of the application.

Order-No.

7071260

7091226

Article

Harvard IonoLine, 2 ml syringe, 3 needle tips
Harvard Needle Tips H20
refill bag with mit 50 needle tips

Harvard Dental International GmbH
Margaretenstr. 2 – 4
15366 Hoppegarten, Germany
Tel: +49 (0) 30/99 28 978-0
Fax: +49 (0) 30/99 28 978-19
info@harvard-dental-international.de
www.harvard-dental.de