Harvard Ionoglas Fill Extra

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
V03.1 03/2018

Harvard Ionoglas Fill Extra

Esthetic self cure glass ionomer cement for restorations

Harvard Ionoglas Fill Extra is a esthetic restorative glass ionomer cement with improved mechanical properties.

Besides it's high fluoride release and high biocompatibility as glass ionomer cement Harvard Ionoglas Fill Extra also has good chemical bonding to dentin and enamel and a tight seal at the dentinal margins. Therefore it requires no enamel and dentin etching. Harvard Ionoglas Fill Extra is radiopaque.

Indications

- Deciduous teeth: final restorative for Class I, II and V (according Black)
- Long term restorative in non-load bearing areas of Class I and II
- Class V restorations
- Intermediate restorations
- Sandwich material for heavy stress bearing Class I and II cavities

Contraindications / side-effects

Do not use Harvard Ionoglas Fill Extra with patients who show an allergy to the material. If allergic reactions occur, further treatment with Harvard Ionoglas Fill Extra must be discontinued and the patient must be advised to consult a physician.

Do not allow the liquid or cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.

Avoid eye contact of the liquid or cement mixture. In case of contact, immediately flush with water and seek medical advise.

Interaction with other materials

Do not mix the powder or liquid of Harvard Ionoglas Fill Extra with any other glass ionomer product.

Application

1. Preparation

Prepare the tooth to be treated in the accustomed manner.

Clean the cavity preparation with pumice and water. Rinse thoroughly and dry, but **do not** desiccate.

If desired, place a matrix band.

By using a cotton pellet apply a dentin conditioner (e.g. Harvard Ionoconditioner) for **20 seconds** to remove the smear layer. Rinse the dentin conditioner with water and dry it in an airstream, but **do not** desiccate.

Pulp capping with Harvard Ionoglas Fill Extra is **contraindicated**. To deep areas close to the pulp apply a small amount of a calcium hydroxide liner (e.g. Harvard CalciumHydroxide).

2. Dosing of powder and liquid

The powder/liquid ratio to achieve a suitable consistency is 3.6 / 1.0. Mix 1 level spoon (orange) of powder and 1 drop of liquid.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

Shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad provided. Avoid compressing powder into the spoon with the inside wall of the bottle.

Turn the liquid bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discount drops** that are obviously not full-sized.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Mixing of powder and liquid

Use a small spatula to rapidly mix all the cement powder into the liquid. Usually one scoop of powder/ one drop of liquid should provide a sufficient amount of mixed cement. The mixed cement should be thixotropic and have a homogenous consistency. Total mixing time is **30 seconds**.

Harvard Ionoglas Fill Extra sticks to metal instruments and should be washed off with cold water prior to setting.

4. Filling

Remove moisture from the cavity surface either with a cotton pellet or very gently with an air syringe blast. It is essential not to desiccate the dentin and enamel surfaces.

Place the mixed cement with a suitable instrument within the working time into the cavity. Please see to it, that no air bubbles will be incorporated.

Form the material with a placement or forming instrument. Optional use a matrix strip to form the surface.

Working time from start of mix at 23°C (74°F): **1:30 minutes**

Net setting time at 37°C (99°F): **3:30 minutes**

Notes

Slightly wet instruments simplify filling and modeling.

Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time reduces adhesion to enamel and dentin.

If a matrix strip is placed, don't remove it before end of net setting time.

After setting, immediately apply **Harvard Ionocoat LC** and light cure for **20 seconds**.

Close bottle immediately after use!

5. Finishing

Final finishing and polishing can begin from about **6 minutes** after start of mixing. Smooth with fine diamond burs, then polish with finishing and polishing discs with graded grain size.

Spray preparation dust away with water and dry the surface in an oilfree airstream.

Apply a thin layer of **Harvard Ionocoat LC** to the final finished surface of the restoration. Immediately light cure for **20 seconds**.

Close bottle immediately after use!

Instruct the patient not to expose the restoration to any pressure for one hour.

Storage information

Storage temperature 4-25°C (39-77°F).

Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Harvard Ionoglas Fill Extra conforms to the DIN EN ISO 9917-1.

Keep the instructions for use for the duration of the application.

Order-No.	Article
	15 g powder / 8 ml liquid, dosage spoon, mixing pad
7052112	Shade A2
7052113	Shade A3
7052135	Shade A3.5

 Manufacturer	Harvard Dental International GmbH Margaretenstr. 2 – 4 15366 Hoppegarten, Germany	 0482 Made in Germany	Tel: +49 (0) 30/99 28 978-0 Fax: +49 (0) 30/99 28 978-19 info@harvard-dental-international.de www.harvard-dental.de	V03.1 03/2018