Harvard Ionoresin Fill LC

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Harvard Ionoresin Fill LC

Light cure glass ionomer cement for restorations

Harvard Ionoresin Fill LC is a light and self cure resin modified glass ionomer restorative cement.

Harvard Ionoresin Fill LC is simple to use. The curing is initiated by common visible light. Light curing allows extended working time and fast setting on demand.

Indications

- Class III and V restorations according to Black, and particularly the restoration of cervical erosions and root surface caries.
- Small class I and II restorations
- Restoration of primary teeth
- Core build-up material
- Intermediate restorations
- Sandwich material for heavy stress bearing Class I and II Cavities

Contraindications / side-effects

- In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used further on.
- Direct and indirect pulp capping

Interaction with other materials

Avoid direct contact with products containing eugenol since eugenol impairs the setting of the material.

1. Tooth Preparation

Clean the cavity preparation with pumice and water. Rinse thoroughly and dry gently, but do not desiccate.

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard CalciumHydroxide) or MTA cement (e.g. Harvard MTA Universal or Harvard MTA XR Flow Fast Opticaps). Refer to the instructions of the product used.

2. Dosing of powder and liquid

The **powder/liquid ratio** to achieve a suitable consistency is **2.8 /1.0**. This can be obtained by mixing 1 level scoop (green) of powder and 2 drops of liquid. A condensable thick mix is required to minimize abrasive loss, but a glossy surface is also need for strong chemical bond strength. Using the appropriate powder / liquid ratio is therefore strongly recommended.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

Shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad provided. **Avoid compressing** powder into the spoon with the inside wall of the bottle.

Turn the liquid bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discount drops** that are obviously not full-sized.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Mixing of powder and liquid

Using a spatula, divide the powder into **2 equal parts**. Spread the liquid across the mixing pad and mix the first half with the whole of the liquid for **15 sec**. Add the second part and mix for **15 sec** to obtain a homogenous mixture. **Total mixing time is 30 sec**.

4. Application of the filling

Mix the required amount of **Harvard Ionoresin Fill LC** and apply it into the prepared cavity by using a suitable placement instrument.

Please see to it, that no air bubbles will be incorporated.

Harvard Ionoresin Fill LC should be placed in the cavity within the **working time (2:00 minutes)** from start of mixing at 23 °C or 74°F). A transparent matrix may be used.

Cure for **20 sec** with a visible light-curing device (wavelength range 400–500 nm, light intensity min. 1000 mW/cm²). Net setting time without any light is approx. 5:30 min.

In case of cavities with a depth of more than 2 mm curing in incremental layers is strongly recommended. Cure each layer separately. In order to improve the self-adhesion, this is very important for the first layer of 2 mm.

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentine.

After Harvard Ionoresin Fill LC has set, remove matrix where applicable and perform finishing under water spray using standard techniques.

Optional a layer of Harvard Ionocoat LC can be applied after the finishing and subsequently light-cured for **20 seconds**.

Warnings

- Unpolymerized material may have an irritant effect and can lead to sensitization against methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, remove the material with absorbent cotton soaked in alcohol and rinse with water. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Harvard Ionoresin Fill LC conforms to the DIN EN ISO 9917-2.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7041216	15 g powder shade A3 / 8 ml liquid dosage spoon, mixing pad

 <p>Manufacturer</p>	<p>Harvard Dental International GmbH Margaretenstr. 2 – 4 15366 Hoppegarten, Germany</p>	 <p>0482 Made in Germany</p>	<p>Tel: +49 (0) 30/99 28 978-0 Fax: +49 (0) 30/99 28 978-19 info@harvard-dental-international.de www.harvard-dental.de</p>	<p>V04 01/2021</p>