Harvard PremiumFill

Harvard Dental International GmbH
 Margaretenstr. 2 – 4
 15366 Hoppegarten, Germany
 Tel: +49 (0) 30/99 28 978-0
 Fax: +49 (0) 30/99 28 978-19
 info@harvard-dental-international.de
 www.harvard-dental.de

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Instructions for Use EN

Universal light cure nano hybrid composite for restorations

Harvard PremiumFill is a light cure nano-filled composite for restorations with high demands for esthetics and durability. Beside its low shrinkage and high wear resistance Harvard PremiumFill shows excellent mechanical properties. And above all Harvard PremiumFill has excellent handling properties for a fast and easy application.

Harvard PremiumFill is based on urethanedimethacrylate resin and inorganic filler particles <1.0 µm. The total filler content is 82 % by weight and 74 % by volume.

Indications

Suitable for all cavity classes.

Contraindications

The placement of Harvard PremiumFill restorations is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard PremiumFill if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Harvard PremiumFill may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calcium hydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol impairs the hardening of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Application

1. Preparation of tooth and color matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Harvard PremiumFill shade. Harvard PremiumFill is shaded according to VITA® Shades.

2. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. Best pulp protection in case of deep excavations is the photo-activated disinfection (PAD). If this is not possible, deep excavations should be covered with a thin layer of calciumhydroxide material.

3. Application of a bonding agent

Condition and apply a bonding agent according to the total etch technique with an etching gel and subsequent application of a light cure bonding agent or with a self etching bonding agent according to the manufacturer instructions. The application of the self etching bonding agent Harvard Interlock (2 components) resp. Harvard Bond SE Mono (1 component) according to the corresponding instructions is recommended.

4. Application and curing of Harvard PremiumFill

Apply Harvard PremiumFill in the selected shade. Place it by using suitable instruments. Harvard PremiumFill is easy to model. Transparent matrix strips may be used. In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement (e.g. Harvard Ionoresin Fill) or a light cure flowable composite (e.g. Harvard PremiumFlow), followed by the application of Harvard PremiumFill. In case of small cavities Harvard PremiumFill can be applied directly.

For an optimum result apply Harvard PremiumFill in layers of max. 2 mm. Light cure each increment separately.

By using a polymerization unit (wavelength range 400 – 500 nm) with a light intensity of at least 1000mW/cm², cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) **20 seconds**
 Darker shades (e.g. A3,5, B3) **30 seconds**
 Hold the light emission window as close as possible to the filling material.

5. Finishing

Contour with fine diamonds, stones or burs. Polish to high gloss with discs or rubber points. Interproximal finishing is accomplished by fine grit finishing strips.

6. Additional informations

- If Harvard PremiumFill in OptiTips® is applied directly in the mouth, for hygienic reasons the OptiTip® is for single use only.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard PremiumFill may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Harvard Dental Internationals sole obligation shall be repair or replacement of the Harvard Dental International product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Only for dental use!

Store product out of reach of children!

Harvard PremiumFill conforms to the ISO 4049.

The instructions for use has to be kept for the duration of the application.

Order-No.	Article
Harvard PremiumFill, 4 g syringe	
7082201	Shade A1
7082202	Shade A2
7082203	Shade A3
7082204	Shade A3,5
7082205	Shade B2
7082206	Shade Den
7082207	Farbe Inc
7082210	Shade Assortment, A1, A2, A3, A3,5, B2, Den, Inc
Harvard PremiumFill, 20 Opti-Tips® ea. 0,3 g	
7083621	Shade A1
7083624	Shade A2
7083625	Shade A3
7083626	Shade A3,5
7083627	Shade B2