**Universelles lichthärtendes
Microhybrid Composite**

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**Universal light cure
microhybrid composite**

Harvard Restore is a light cure microhybrid composite for anterior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05-1.5 µm. The total filler content is 81 % by weight and 65 % by volume.

Indications
Suitable for all cavity classes.

Contraindications
The placement of Harvard Restore restorations is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard Restore if the patient is known to be allergic to any of the ingredients.

Side effects
In singular cases, Harvard Restore may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

Incompatibility with other materials
Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor allow not polymerized composite to come into contact with materials containing eugenol.

Application
1. Preparation of tooth and color matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching. Ascertain the tooth shade while teeth are still moist and select the appropriate Harvard Restore shade.

Harvard Restore is shaded according to VITA®-Shades.

2. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. For pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material (e.g. Harvard CalciumHydroxide).

3. Application of a bonding agent

Condition and apply a bonding agent according to the total etch technique (e.g. Harvard Etch) and subsequent application of a light cure bonding agent (e.g. Harvard Bond TE Mono).

Alternatively a self etching bonding agent e.g. Harvard InterLock (two component adhesive) resp. Harvard Bond SE Mono (one component adhesive) can be applied according to the corresponding instructions.

4. Application and curing of Harvard Restore

Apply Harvard Restore in the selected shade. Place it by using suitable instruments. Harvard Restore is easy to model. Matrix strips may be used. In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement (e.g. Harvard Ionoglas Fill LC) or a light cure flowable composite (e.g. Harvard RestoreFlow), followed by the application of Harvard Restore.

In case of small cavities Harvard Restore can be applied directly.

For an optimum result apply Harvard Restore in layers of max. 2 mm. Light cure each increment separately.

By using a polymerization unit (wavelength range 400 – 500 nm) with a light intensity of at least 1000 mW/cm², cure each increment as follows:
Lighter shades (e.g. A1, A2, B2) **20 seconds**
Darker shades (e.g. A3.5, B3) **30 seconds**
Hold the light emission window as close as possible to the filling material.

5. Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs. Polish to high gloss e.g. with polishing discs or polishing brushes.

6. Additional informations

- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard Restore may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!
Store product out of reach of children!

Harvard Restore conforms to the DIN EN ISO 4049. Keep the instructions for use for the duration of the application.

Order-No.	Article
Harvard Restore, 4 g syringe	
7083201	Shade A1
7083202	Shade A2
7083203	Shade A3
7083204	Shade A3,5
7083221	Shade B1
7083222	Shade B2
7083215	Shade Inc
7083209	Shade A1 opaque
7083210	Shade A2 opaque

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