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Instructions for Use EN

Vinylpolysiloxane impression material, heavy bodied

Harvard SuperSnap® Heavy Regular Snap is a super hydrophilic heavy bodied impression material:

- for the double-mix technique (1 step/2 phase) with Harvard SuperSnap® Light Regular Snap or Harvard SuperSnap® Medium Regular Snap
- for the correction impression technique (2 step/2 phase) with Harvard SuperSnap® Light or Harvard SuperSnap® Medium

It is characterized by a long working time in conjunction with a short intraoral setting time (super snapset).

Harvard SuperSnap® Heavy is a precision impression material based on addition curing vinyl silicones delivered in automix-cartridges. It can be extruded directly and bubble-free onto the tray.

Harvard SuperSnap® Heavy shows strong hydrophilic properties. Therefore it improves the wettability and yields an exact reproduction of subgingival areas and other moist surfaces.

Harvard SuperSnap® Heavy meets the requirements of ISO 4823, type 2 (medium-bodied consistency).

Contraindications/side effects

None known.

Application

1. Preparation of the automix-cartridge

Harvard Dental International GmbH only recommends for Harvard SuperSnap® Heavy the use of mixing tips type Harvard Auto-T 1:1 O-green.

Mount the cartridge into the application gun. Remove and dispose the closure-cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved. Install the mixing tip onto the cartridge and turn it 90° clockwise until it locks in position and extrude needed material.

Note:

Store used cartridge with fixed used mixing tip. **Do not re-use the cartridge cap.** Closing the cartridge by turning on the cap between use could cause cross contamination of components and thus premature setting. Leave the used mixing tip on the cartridge. It serves as a cap.

2. Impression taking

Double-mix technique

Load the impression tray with Harvard SuperSnap® Heavy Regular Snap and apply Harvard SuperSnap® Light Regular Snap or Harvard SuperSnap® Medium Regular Snap around the preparation. Tray loading and syringing around the preparation have to be completed at the same time.

Insert the tray not later than **1:45 minutes** from start of mixing (working time at 23°C or 74°F) into the mouth. Maintain position for at least **1:45 minutes**.

Two-step (correction) impression technique

Load the impression tray with Harvard SuperSnap® Heavy Regular Snap. Insert the loaded tray in the mouth not later than **1:45 minutes** (working time at 23°C or 74°F) after start of mix. Maintain position for at least **1:45 minutes**.

Cut drain channels into the pre-impression.

Load the auto mixed Harvard SuperSnap® Light or Harvard SuperSnap® Medium onto the tray with the preliminary impression and/or directly around the teeth. Position the tray within the working time of the correction impression material into the mouth. Maintain position until hardening has finished. Refer to the corresponding user instructions.

3. Casting

The impression should be stored under normal room conditions (max. 25°C or 77°F). The impression can be poured 30 minutes after removal from the mouth. No other time limitation does exist.

4. Electroplating

All impressions can be copper- or silverplated.

5. Disinfection

Impressions can be disinfected with standard disinfection solutions.

Notes:

Do not contact the product with latex gloves, they may inhibit the setting. Do not touch teeth with latex gloves.

Do not use with standard (condensation reaction type) silicone impression materials.

Storage

Do not store above 25 °C (77 °F). Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children! Keep the instructions for use for the duration of the application.

Order-No.	Article
7083831	Harvard SuperSnap® Heavy, Regular Snap, 2 automix cartridges ea. 50 ml, 12 mixing tips
7098010	Harvard Auto-T 1:1 O-green, refill bag with 50 mixing tips
7095100	Harvard Dispenser Automix 1:1

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