# Harvard SuperSnap® Medium Quick Snap

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## Instructions for Use EN

### Vinylpolysiloxane impression material, light bodied

Harvard SuperSnap® Medium Quick Snap is a super hydrophilic low viscosity impression material:

- for the double-mix technique (1 step/2 phase) with Harvard SuperSnap® Putty Soft Quick Snap or Harvard SuperSnap® Heavy Quick Snap
- for the correction impression technique (2 step/2 phase) with Harvard SuperSnap® Putty
- for relining, functional impressions, impressions for denture repairs and for impressions using the copper band technique

It is characterized by a long working time in conjunction with a short intraoral setting time (super snapset).

Harvard SuperSnap® Medium is a precision impression material based on addition curing vinyl silicones delivered in automix-cartridges. It can be extruded directly onto the teeth or the tray.

Harvard SuperSnap® Medium is thixotropic. That means it stays in position in the tray, but under pressure it flows. The product shows strong hydrophilic properties. Therefore it improves the wettability and yields an exact reproduction of subgingival areas and other moist surfaces.

Harvard SuperSnap® Medium meets the requirements of ISO 4823, type 3 (light-bodied consistency).

#### Contraindications/side effects

None known.

### Application

#### 1. Preparation of the automix-cartridge

Harvard Dental International GmbH only recommends for Harvard SuperSnap® Medium the use of mixing tips type Harvard Auto-T 1:1 O-yellow.

Mount the cartridge into the application gun. Remove and dispose the closure-cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved. Install the mixing tip onto the cartridge and turn it 90° clockwise until it locks in position and extrude needed material. For direct extruding around the teeth use Intra Tips yellow.

**Note:**  
Store used cartridge with fixed used mixing tip. **Do not re-use the cartridge cap.** Closing the cartridge by turning on the cap between use could cause cross contamination of components and thus premature setting. Leave the used mixing tip on the cartridge. It serves as a cap.

### 2. Impression taking Double-mix technique

Load the impression tray with the tray impression material (with Harvard SuperSnap® Putty Soft Quick Snap or Harvard SuperSnap® Heavy Quick Snap) and apply Harvard SuperSnap® Medium Quick Snap around the preparation. Tray loading and syringing around the preparation have to be completed at the same time. Insert the tray not later than **1:10 minutes** from start of mixing (working time at 23°C or 74°F) into the mouth. Maintain position for at least **1:15 minutes**.

### Two-step (correction) impression technique

Cut drain channels into the pre-impression. Load the auto mixed Harvard SuperSnap® Medium Quick Snap onto the tray with the preliminary impression and/or directly around the teeth. Insert the tray not later than **1:10 minutes** from start of mixing (working time at 23°C or 74°F) into the mouth. Maintain position for at least **1:15 minute**.

### 3. Casting

The impression should be stored under normal room conditions (max. 25°C or 77°F). The impression can be poured 30 minutes after removal from the mouth. No other time limitation does exist.

### 4. Electroplating

All impressions can be copper- or silverplated.

### 5. Disinfection

Impressions can be disinfected with standard disinfection solutions.

#### Notes:

Do not contact the product with latex gloves, they may inhibit the setting. Do not touch teeth with latex gloves.

Do not use with standard (condensation reaction type) silicone impression materials.

#### Storage

Do not store above 25 °C (77 °F). Do not use after expiry date.

#### Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

#### Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory

asserted, including warranty, contract, negligence or strict liability.

#### For dental use only!

**Store product out of reach of children!**  
Keep the instructions for use for the duration of the application.

| Order-No. | Article                                                                                                        |
|-----------|----------------------------------------------------------------------------------------------------------------|
| 7083812   | Harvard SuperSnap® Medium, Quick Snap, green, 2 automix cartridges ea. 50 ml, 12 mixing tips 1:1, 4 Intra Tips |
| 7098000   | Harvard Auto-T 1:1 O-yellow, refill bag with 50 mixing tips for Light/Medium – Regular/Quick                   |
| 7083619   | Harvard Intra Tips yellow refill bag with 50 intra tips.                                                       |
| 7095100   | Harvard Dispenser Automix 1:1                                                                                  |

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