Harvard BioPOWER Cem HandMix

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Harvard BioPOWER Cem HandMix

Bioactive dual cure MTA-modified composite-hybrid cement for luting of crowns and bridges

Product description

Harvard BioPOWER Cem is a dual cure resin modified MTA glass ionomer luting cement.

Indications/Intended use

- Luting of crowns and bridges (e.g. zirconia, metal, ceramic)
- Luting of posts

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

Do not use this product in cases where known allergies to the components exist.

Harvard BioPOWER Cem is not indicated for luting of inlays, onlays, crowns and bridges made of composites.

Due to its high opacity, Harvard BioPOWER Cem should not be used if the restoration has to meet high esthetic requirements.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Interaction with other materials

Avoid direct contact with products containing eugenol since eugenol impairs the setting of the material.

Application

1. Tooth Preparation

Clean the luting side areas of the restoration according to manufacturer instructions.

2. Dispensing and Mixing

The **powder/liquid ratio** is. **4.0/1.0**. This can be obtained by mixing **1 level of the supplied scoop of powder and 2 drops of liquid**.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

For accurate dispensing of Harvard BioPOWER Cem powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon with the inside wall of the bottle.

For dispensing of Harvard BioPOWER Cem liquid turn the bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discard drops** that contain bubbles and are obviously not full-sized.

Use a small spatula to rapidly mix all the cement powder into the liquid.. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **60 seconds**.

Harvard BioPOWER Cem is extremely moisture sensitive!
Immediately after use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Cementing

Mix the required amount of Harvard BioPOWER Cem and apply it with a suitable placement instrument into the lumens of the crown or bridge (approx. half) resp. to the bonding surface of inlays or onlays in a thin layer. Seat the restoration immediately (**working time 2:00 minutes** from start of mixing at 23°C). Net setting time without any light is **3:00 minutes**.

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified (**3:00 minutes**). Setting can be optimized with **20 seconds** light cure with a suitable dental light cure unit (wavelength range 400–500 nm, light intensity min. 1000 mW/cm²).

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentine.

4. Additional Notes/Warnings

- Unpolymerized material may have an irritant effect and can lead to sensitization against methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

5. Storage

Store at 10 - 25°C (50 - 77°F). Store protected from moisture. After first opening use up within 3 months and before the end of expiry date. Do not use after expiry date.

Composition

Powder: MTA, Zincoxide, dental glass, photo initiators
Liquid: Polyacrylic acid, methacrylates

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order-No. Article

7071421 5 g powder / 3 ml liquid, dosage spoon, mixing pad universal