Harvard POWERmer Bulk HandMix

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Harvard POWERmer Bulk HandMix

High strength, self adhesive, dual cure resin modified glass ionomer cement for restorations

Product description

Harvard POWERmer Bulk is a dual cure resin reinforced self-adhesive glass ionomer restorative material with high strength that is also suitable as an amalgam alternative. Harvard POWERmer Bulk has good bonding properties to dentine and enamel and does not shrink. With Harvard POWERmer Bulk, the cavity can be filled in just one layer up to the occlusal surface (Bulk-Fill). No adhesive is required.

Indications/Intended use

- Restorations of cavity classes I, II, III and V
- Core build-up
- Base and liner

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

Do not use this product in cases where known allergies to the components exist. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Interaction with other materials

Avoid direct contact with products containing eugenol since eugenol impairs the setting of the material.

Application

1. Tooth Preparation

Prepare the cavity. Rinse thoroughly with water and dry gently, but do not desiccate.

For pulp protection areas close to the pulp should be covered with a thin layer of an alkaline liner (e.g. Harvard BioCal®-Line).

2. Dispensing and Mixing

The **powder/liquid ratio** is **4.3/1.0**. This can be obtained by mixing **1 level of the supplied scoop of powder** and **1 drop of liquid**.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

For accurate dispensing of Harvard POWERmer Bulk powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon with the inside wall of the bottle.

For dispensing of Harvard POWERmer Bulk liquid turn the bottle vertically with the tip about 5cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discard drops** that contain bubbles and are obviously not full-sized.

Use a small spatula to rapidly mix all the cement powder into the liquid.. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **60 seconds**.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

The powder bottle contains a dessicant pellet. Do not damage it!

3. Filling and Finishing

Mix the required amount of Harvard POWERmer Bulk and apply it into the prepared cavity by using a suitable placement instrument. Fill in a single layer up to the desired height.

Please see to it, that no air bubbles will be incorporated!

Fill the material into the cavity within **1:30 minute** (working time from start of mixing at 23 °C or 74°F) and shape the surface with suitable instruments. A matrix can be used to shape the surface.

Cure **for 20 seconds** with a visible light-curing device (wavelength range 400–500 nm, light intensity min. 1000 mW/cm²). Net setting time without any light is **3:30 – 4:00 minutes**.

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the dental enamel and the dentine.

After Harvard POWERmer Bulk has set, remove matrix where applicable and perform finishing under water spray using standard techniques and polish.

4. Additional Notes/Warnings

- Unpolymerized material may have an irritant effect and can lead to sensitization against methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

5. Storage

Store at 10 - 25°C (50 - 77°F). Store protected from moisture. Do not use after expiry date.

Composition

Powder: Dental glass, photo initiators

Liquid: polyalkenoate acid, methacrylates

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbHs sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order-No.

7071222

Article

15 g powder / 7 ml liquid, dosage spoon, mixing pad A2

 Manufacturer	Harvard Dental International GmbH Robert-Bosch-Strasse 2 25335 Elmshorn, Germany	CE 0482 MD Made in Germany	Tel: +49 (0) 30/99 28 978-0 Fax: +49 (0) 30/99 28 978-19 info@harvard-dental-international.de www.harvard-dental.de	V02_MDR 2025-04
---	--	--	--	-------------------------------------