## Harvard VeneerCem

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
## Harvard VeneerCem

### Shade-adaptive light cure veneer composite cement for translucent veneers

#### Product description

Harvard VeneerCem is a light cure flowable composite cement for luting of translucent veneers.

#### Indications/Intended use

Permanent cementation of translucent veneers made from ceramic and composite

#### Performance features

The performance features of the product meet the requirements of the intended use.

#### Contraindications

**In case of known hypersensitivity to methacrylates this product should not be used.**

Do not use Harvard VeneerCem for teeth with severe discoloration.

#### Patient target group

Persons who are treated during a dental procedure.

#### Intended users

This medical device should only be used by a professionally trained dental practitioner.

#### Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

#### Application

**Before starting work, the luting side areas of the veneers must already be prepared according to the manufacturer's instructions.**

##### 1. Preparation

Remove possible temporary restoration and clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water). Rinse with water and dry with oil-free air. Avoid overdrying.

##### 3. Preparation of the tooth

Isolate the tooth to be veneered with interproximal strips to protect adjacent teeth (not being veneered) from the etchant and bonding agent.

Apply an etchant (e.g. Harvard Etch) onto the enamel areas to be veneered. Leave the etchant in place for 15 - 30 seconds, and rinse with plenty of water. After drying place a thin coat of suitable bonding agent as described in the corresponding user instructions.

##### 4. Application of Harvard VeneerCem

Screw off the cap of the syringe and substitute it by one of the supplied NeedleTips.

After use discard the NeedleTip and attach the cap to the syringe again. NeedleTips are for single use only due to hygienic reasons.

Apply Harvard VeneerCem in an even thin layer to the luting side areas of the veneers.

Remove matrix strips prior to placing veneers on teeth. This step assures complete and passive seating of the veneer, even in multiples.

**Seat the veneer and hold it in place during excess removal.**

Immediately after placing the veneer, remove all excess luting material with a microbrush/brush/dental floss or another suitable instrument. Make sure to remove excess material especially in areas that are difficult to reach (proximal areas, gingival margins).

#### Note:

To avoid building of an inhibition layer during polymerization cover the

margins of the veneers with a glycerin gel.

Light cure Harvard VeneerCem in segments, starting with the proximal margins.

By using a polymerization unit (wavelength range 400-500 nm) with a light intensity of at least 1000mW/cm<sup>2</sup>, cure each segment for **20 seconds** through veneer material and tooth structure.

Rinse off glycerin gel.

After polymerization finish and polish the margins as usual.

#### Storage

Do not store above 25°C (77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

#### Additional Notes/Warnings

- The ambient light of the dental lamp may start polymerization of the composite.
- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized composite may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

#### Composition

Dimethacrylates, prepolymerized resin, dental glass, silicon dioxide, photo initiators.

#### Disposal

Disposal of the product according to local authority regulations.

#### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

#### Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed/>).

#### Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

#### Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

#### Caution:

**Federal law restricts the sale of this device to or by the order of a dentist.**

Keep the instructions for use for the duration of the application.

#### Order-No. Article

7100020	Harvard VeneerCem 2 x 1 ml syringe, 8 needle tips H18
7095158	Harvard NeedleTips H18 refill bag with 50 needle tips.