Harvard Bond SE Dual

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Harvard Bond SE Dual

Self-etching, dual-curing adhesive (two component system)

Product description

Harvard Bond SE Dual is a simple to use **self-etching dual cure adhesive** for a durable and strong bonding of **composites, compomers and resin modified glass ionomer cements to enamel and dentin**. Harvard Bond SE Dual works optimally under slightly wet conditions. Harvard Bond SE Dual consists of the components PART A and PART B that were mixed before application. It is bonding to **dual cure, self cure and light cure** composite restorative materials.

Indications/Intended use

Adhesive for:

- Core build-ups and cementing of posts with dual and self cure composites (e.g. Harvard ZirconCore)
- Cementing of posts with dual and self cure composite cements
- Cementing of inlays, onlays, crowns and bridges with dual and self cure composite cements
- Restorations with dual cure, self cure and light cure composites

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used.

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of Harvard Bond SE Dual. Neither store the material in proximity of eugenol containing products, nor let the material allow coming into contact with materials containing eugenol.

Application

1. Isolation

Rubber dam is the recommended method of isolation.

2. Cavity Preparation

Prepare the cavity according to the principles of the adhesive restoration technique.

3. Pulp Protection

For pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material material (e.g. Harvard BioCal®-CAP or Harvard CalciumHydroxide).

4. Application of Harvard Bond SE Dual

One drop of Harvard Bond SE Dual PART A and one drop of Harvard Bond SE Dual PART B were combined in a mixing pallet and mixed for 5-10 seconds.

Notes:

Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.

4.1. Application for Core Build-ups, Cementing of Inlays, Onlays, Crowns and Bridges and for Composite based Restorations

Apply the homogeneous mixture generously with a brush onto the slightly wet enamel- and dentin surfaces (as they are usual in clinical procedures) for **30 seconds** with agitation. The material should build a homogeneous layer. Air thin for **10 seconds** to remove the volatile components and to disperse the adhesive. Then light cure for **20 seconds** with a dental halogen light unit or an LED (wavelength 400–500 nm, light intensity min. 1000 mW/cm²) and place the restorative material.

For a maximal adhesion it is strongly recommended to use the light cure mode.

If light cure is absolutely impossible, the adhesive will also cure in the autocure mode. After application of the adhesive air thin to remove all volatile components. There must remain a sticky layer. Then apply the dual cure, self cure or light cure restorative material directly.

4.2. Application in the root canal (cementing of posts)

Prepare and clean the root canal with e.g. sodium hypochlorite solution, rinse and remove excess solution from the canal with a soft paper point.

Apply the homogeneous mixture of Harvard Bond SE Dual Part A and Harvard Bond SE Dual Part B generously with a brush onto the slightly wet root canal walls for 15 seconds with agitation. Repeat procedure 1 – 2 times. All the dentin surfaces must kept wet with the primer over the specified time.

Remove excess Harvard Bond SE Dual with dry paper points. Dry the root canal cautiously with oilfree air for 15 seconds to remove all volatile components and to disperse the adhesive to an even layer. Light cure all areas that are available for a dental curing unit for 20 seconds.

5. Restorative Placement, Cure and Finishing

Refer to manufacturer instructions for placement of the dual cure, self cure or light cure restorative material.

Storage

Do not store above 25 °C (77 °F)! Protect from direct sunlight. Do not use after expiry date.

Warnings

- Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.
- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Composition

Hydrophilic methacrylates, MDP, photo initiators, catalysts, photo initiators.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order-No.

7083602

7083605

7095156

Article

Harvard Bond SE Dual 2x5 ml bottle, 50 Microbrush®, 5 mixing wells

Harvard Bond SE Dual Refill, 2x5 ml bottle

Microbrush®, refill bag with 50 applicators