Harvard InterLock® ONE Self Cure Activator

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Instructions for Use EN

Self Cure Activator for use in combination with Harvard InterLock® ONE Universal Adhesive

Product description

Harvard InterLock® ONE Self Cure Activator is designed as a product for additional and/or extended self etch applications of Harvard InterLock® ONE. It is applied in combination with Harvard InterLock® ONE in all the special cases, where light cure of the adhesive cannot be ensured (e.g. in the root canal).

The mixture of Harvard InterLock® ONE/ Harvard InterLock® ONE Self Cure Activator is suitable for moist, wet and dry surfaces.

Indications/Intended use

Activator for Harvard InterLock® ONE for the self-etch technique:

- Adhesive cementing of posts (and possibly subsequent core build-up) with dual or self cure composites (e.g. Harvard ZirconCore).

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

In singular cases, Harvard InterLock® ONE/Harvard InterLock® ONE Self Cure Activator may cause a sensitizing. In these cases, the material should not be further used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Incompatibility with Other Materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Application

1. Isolation

Rubber dam is the recommended method of isolation.

2. Preparation of the Root Canal

Prepare the root canal appropriate to the selected post.

3. Application of the Adhesive

Mix one drop of Harvard InterLock® ONE and one drop of Harvard InterLock® ONE Self Cure in a mixing pallet for **5-10 seconds** under subdued light.

Notes:

Do not interchange lids of the bottles, because this can lead to a cross-contamination of the liquids.

Apply the Harvard InterLock® ONE/Harvard InterLock® ONE Self Cure Activator mixture with a suitable brush onto the root canal walls for **30 seconds** with agitation. Applying the adhesive on the surface without agitation is inadequate. Repeat procedure 1 – 2 times. The material should build a homogeneous layer.

Remove excess material **thoroughly** e.g. with paper points. Air thin gently (**at least 5-10 seconds**) to remove the volatile components and to disperse the adhesive until a glossy and uniform layer results.

For optimal adhesion it is recommended to light cure the coronal part of the root canal with a suitable dental halogen light unit or an LED (wavelength 400–500 nm, light intensity min. 1000 mW/cm²) for **20 seconds**.

4. Post Cementing

Prepare the post according to manufacturer instructions.

Place the post with a flowable dual cure or self cure composite (e.g. Harvard ZirconCore) into the root canal. Refer for this (and for core build-up) to manufacturer instructions of the dual or self cure composite.

Additional informations/Warnings

Unpolymerized material may have an irritating effect and may lead to a sensitizing reaction against methacrylates.

- Avoid contact with skin, mucous membrane and eyes.
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Composition

Catalysts, photoinitiators.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct,

indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7083611	Harvard InterLock® ONE 5 ml bottle
7083613	50 SingleDose ea. 0,1 ml
7083612	Self Cure Activator 5 ml bottle

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