## Harvard IonoResin Cem Extra

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## Harvard IonoResin Cem Extra

### Light and self cure resin-modified glass ionomer luting cement

#### Product description

**Harvard IonoResin Cem Extra** is a light cured resin reinforced glass ionomer luting and lining cement.

Beside its excellent fluoride release and the good mechanical properties Harvard IonoResin Cem Extra has also a very low solubility and a tight seal at the dental margins. Because of its radiopacity it ensures easy postoperative diagnosis.

**Harvard IonoResin Cem Extra OptiCaps®** can quickly and simply be activated and extruded with the Harvard Applier OptiCaps®. Capsule mixing (mixing time 10 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/min. Application can be done directly from the capsule.

#### Indications/Intended use

Cementing of crowns, bridges, inlays and onlays made from metal, resin, ceramic fused to metal or ceramic)

#### Performance features

The performance features of the product meet the requirements of the intended use.

#### Contraindications

- In singular cases, the material may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used further on.
- Direct and indirect pulp capping

#### Patient target group

Persons who are treated during a dental procedure.

#### Intended users

This medical device should only be used by a professionally trained dental practitioner.

#### Interaction with other materials

Avoid direct contact with products containing eugenol since eugenol impairs the setting of the material.

Do not mix the powder or liquid of Harvard IonoResin Cem Extra with any other glass ionomer product.

#### Application

##### 1. Tooth preparation

Prepare the tooth to be restored in the usual manner.

For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard BioCal® CAP, Harvard CalciumHydroxide) or MTA cement (e.g. Harvard MTA Universal or Harvard MTA XR Flow Fast Opticaps). Refer to the corresponding instructions for use.

##### 2. Dosing of powder and liquid (only for HandMix)

The **powder/liquid ratio** to achieve a suitable consistency is **2.0 /1.0**. This can be obtained by mixing **1 level (orange) scoop of powder and 3 drops of liquid**.

#### Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

For accurate dispensing of Harvard IonoResin Cem Extra powder shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder by using the scraper at the top of the bottle and carry it onto the mixing pad. Avoid compressing powder into the spoon.

If any bubbles are present in the Harvard IonoResin Cem Extra liquid, lightly tap the bottle with the fingers holding it. Discard drops that contain bubbles and are obviously not full-sized.

For dispensing of Harvard IonoResin Cem Extra liquid turn the bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

#### 3. Mixing

##### 3.1. Mixing of powder and liquid (HandMix)

Use a small spatula to rapidly mix all the cement powder into the liquid. The mixed cement should be thixotropic and have a smooth consistency and glossy appearance. Total mixing time is **30 seconds**.

##### 3.2. Mixing Harvard IonoResin Cem Extra OptiCaps®

See "Instructions for the activation and mixing of Harvard IonoResin Cem Extra OptiCaps®"

**Mixing time for the OptiCaps® is 10 seconds.**

Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting which may impair or prevent the application of the material.

All the cement required should be dispensed from the OptiCaps® within **30 seconds** from the start of mixing.

#### 4. Cementing

Avoid water and saliva contamination during application and setting of the cement. For a dry working area adequate isolation of the tooth is required.

Prepare the restoration according to the manufacturer instructions.

Apply the mixed cement into the lumens of the crown or bridge (approx. half) resp. to the bonding surface of inlays or onlays in a thin layer. Insert the restoration immediately.

**Working time** from start of mix at 23 °C (74 °F)

**Handmix: 2:30 minutes**

**OptiCaps®: 2:00 minutes**

Net setting time without any light is approx. 4:00 minutes (at 37°C / 99 °F).

Remove excess cement at the first setting stage.

Maintain isolation until the set of the cement is verified (ca. 4:00 minutes). Setting can be optimized with **20 seconds** light cure with a suitable dental light cure unit (wavelength range 400–500 nm, light intensity min. 1000 mW/cm<sup>2</sup>).

Note: Higher temperatures will shorten the working time, lower temperatures will prolong the working time.

An overextended working time will cause the loss of adhesion to the enamel and the dentine.

#### Notes/Warnings

- The material can also be used as liner under composite restorations.
- Unpolymerized material may have an irritant effect and can lead to sensitization against methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, remove the material with absorbent cotton soaked in alcohol and rinse with water. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

#### Composition

Powder: Dental glass, photo initiators

Liquid: polyalkenoate acid, methacrylates

#### Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date.

Harvard IonoResin Cem Extra OptiCaps® are for single use only.

#### Disposal

Disposal of the product according to local authority regulations.

#### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

#### Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

#### Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

#### Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Harvard IonoResin Cem Extra conforms to the DIN EN ISO 9917-2.

Keep the instructions for use for the duration of the application.

#### Order-No.

7061116

7061251

7092000


#### Article

15 g powder, shade Universal / 10 ml liquid, dosage spoon, mixing pad


50 OptiCaps® ea. 0.4 g, shade Universal

Harvard Applier OptiCaps®


 Manufacturer	<b>Harvard Dental International GmbH</b> <b>Robert-Bosch-Strasse 2</b> <b>25335 Elmshorn, Germany</b>	 0482  Made in Germany	Tel: +49 (0) 30/99 28 978-0 Fax: +49 (0) 30/99 28 978-19 <a href="mailto:info@harvard-dental-international.de">info@harvard-dental-international.de</a> <a href="http://www.harvard-dental.de">www.harvard-dental.de</a>	V02_MDR 2025-10
---	---	---	---	--------------------

**Instructions for the activation and mixing of Harvard IonoResin Cem Extra OptiCaps®**

Harvard OptiCaps® before activation.


For activation of the Harvard OptiCaps®, press the plunger on a hard and plane surface to the end into the capsule.


Insert the Harvard OptiCaps® into the Harvard Applier OptiCaps® and click once.


Note: The plunger must be at the same level as the bottom of the capsule.


Insert the Harvard OptiCaps® into a mixer, close lid and mix immediately for **10 seconds** at 4,300 oscillations/minute.

Insert the Harvard OptiCaps® into the Harvard Applier OptiCaps®. Remove the pin from the nozzle. If you leave it, the capsule can burst.

Pull the lever 2 times (2 clicks) to prime the the Harvard OptiCaps®.

Extrude the mixed material on a glass plate or apply directly. Unlock the Harvard Applier OptiCaps® (press button A) and remove the Harvard OptiCaps®. The optimal amount of mixed material is guaranteed only with the Harvard Applier OptiCaps® (Order-No. 7092000).

For the selection of a suitable capsule mixer, our sales and marketing colleagues are gladly available to you.