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Dosierlöffel, Mischblock

Harvard Ionoglas Fill

Classical Glass Ionomer Restorative Cement

Product description

Harvard Ionoglas Fill is a non sticky condensable esthetic restorative glass ionomer cement. Besides its high fluoride content and excellent biocompatibility Harvard Ionoglas Fill has also good chemical bonding to dentin and enamel. Therefore it requires no enamel and dentin etching. The early resistance to water uptake permits one visit treatment. Because of its radiopacity it ensures easy postoperative diagnosis.

Indications/Intended use

- Deciduous teeth: final restorative for Class I, II and V (according Black)
- Long term restorative in non-load bearing areas of Class I and II
- Class V restorations
- Intermediate restorations
- Sandwich material for heavy stress bearing Class I and II cavities
- Core build-up material

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

- Pulp capping
- Allergic reactions to glass ionomer cements

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Application

1. Preparation

Prepare the tooth to be treated in the accustomed manner.

Clean the cavity preparation with pumice and water. Rinse thoroughly and dry, but **do not** desiccate.

If desired, place a matrix band.

By using a cotton pellet apply a dentin conditioner to remove the smear layer. Refer to the instructions of the product used. Rinse the dentin conditioner with water and dry it in an airstream, but **do not** desiccate.

For pulp protection areas close to the pulp should be covered with a thin layer of pulp capping material (e.g. Harvard BioCal®-CAP or Harvard CalciumHydroxide).

2. Dosing of powder and liquid

The powder/liquid ratio to achieve a suitable consistency is 3.6 / 1.0. Mix 1 level spoon (orange) of powder and 1 drop of liquid.

Note:

The spoon dosage is only an approximate reproduction of the nominal mixing ratio.

A condensable thick mix is required to minimize abrasive loss, but a glossy surface is also needed for strong chemical bond strength. Using the appropriate powder / liquid ratio is therefore strongly recommended.

Shake the bottle to loosen the powder. Overfill the spoon with the powder, level the powder for using the collar at the top of the bottle and carry it onto the mixing pad provided. Avoid compressing powder into the spoon with the inside wall of the bottle.

Turn the liquid bottle vertically with the tip about 5 cm above the mixing pad. Steady your hand and squeeze the bottle gently to dispense one drop. If any bubbles are present, lightly tap the bottle with the fingers holding it. **Discard drops** that are obviously not full-sized.

After use, tightly close both liquid and powder bottles to prevent exposure to moisture.

3. Mixing of powder and liquid

Using a spatula, divide the powder into 2 equal parts. Spread the liquid across the mixing pad and mix the first half with the whole of the liquid for **15 seconds**. Add the second part and mix for **15 seconds** to obtain a homogenous mixture. **Total mixing time is 30 seconds**.

4. Filling

Remove moisture from the cavity surface either with a cotton pellet or very gently with an air syringe blast. It is essential not to desiccate the dentin and enamel surfaces. Place the mixed cement with a suitable instrument within the working time into the cavity. Please see to it, that no air bubbles will be incorporated.

Form the material with a placement or forming instrument. Optional use a matrix strip to form the surface.

Working time from start of mix at 23°C (74°F): **1:20 minutes**
Net setting time at 37°C (99°F): **4:00 minutes**

Notes

Slightly wet instruments simplify filling and modeling. Higher temperatures will shorten the working time, lower temperatures will prolong the working time. An overextended working time reduces adhesion to enamel and dentin.

If a matrix strip is placed, don't remove it before end of net setting time.

Immediately after setting a suitable varnish can be applied.

5. Finishing

Final finishing and polishing can begin from about **6 minutes** after start of mixing. Smooth with fine diamond burs, then polish with finishing and polishing discs with graded grain size. Spray preparation dust away with water and dry the surface in an oilfree airstream.

If desired apply a suitable varnish to the final finished surface of the restoration.

Instruct the patient not to expose the restoration to any pressure for one hour.

6. Additional Notes/Warnings

- Do not use the material with patients who show an allergy to the material. In case of allergic reactions immediately stop the application, and advise the patient to consult a physician. An operator, who has a history of allergy to glass ionomer cements, should not handle with this material.
- Do not allow the cement mixture to contact the oral tissues or skin. In case of contact, remove the material with absorbent cotton soaked in alcohol and rinse with water.
- Avoid eye contact of the cement mixture. In case of contact, immediately flush with water and seek medical treatment.
- Do not mix the powder or liquid of Harvard Ionoglas Fill with any other glass ionomer product.
- Keep away from children!

Storage

Storage temperature 4-25°C (39-77°F). Store protected from moisture. Do not use after expiry date.

Composition

Dental glass, polyalkenoate acids.

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Harvard Ionoglas Fill conforms to the DIN EN ISO 9917-1.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7051110	10 g powder, shade A2 / 5,6 ml liquid, dosage spoon, mixing pad
7051115	15 g powder, shade A2 / 8 ml liquid dosage spoon, mixing pad
7051116	15 g powder, shade A3 / 8 ml liquid dosage spoon, mixing pad
7051120	35 g powder, shade A2 / 20 ml liquid dosage spoon, mixing pad