Harvard

MTA-CAP OptiCaps®

(MTA XR Flow Fast)

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Instructions for Use EN

Bioactive MTA cement for pulp capping, flowable, extra radiopacity, extra fast

Product description
Harvard MTA-CAP OptiCaps® is a biocompatible endodontic cement in OptiCaps®. Due to it's unique flowable consistency in combination with a short setting time Harvard MTA-CAP is the ideal choice for **pulp capping**. Harvard MTA-CAP powder consists of very fine hydrophilic particles of several hydraulic mineral oxides. After contact with the liquid it forms a gel that hardens to an impermeable barrier. Harvard MTA-CAP OptiCaps® are easily activated and the content of the capsule is easily ejected with the Harvard Applier OptiCaps®. Mixing of Harvard MTA-CAP OptiCaps® (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/minute.

Indications/Intended use
• Pulp capping (direct and indirect)

Performance features
The performance features of the product meet the requirements of the intended use.

Contraindications
None known

Patient target group
Persons who are treated during a dental procedure.

Intended users
This medical device should only be used by a professionally trained dental practitioner.

Application

1. Activation and mixing
See "Instructions for the activation and mixing of Harvard OptiCaps®"
Mixing time for the Harvard MTA-CAP OptiCaps® is 30 seconds.
Attention:
Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting and lag times may impair or prevent application of the material. The material must be extruded within **10 seconds** after the end of mixing. To prevent dehydration during setting, apply Harvard MTA-CAP immediately after mixing.

2. Working time
Working time at 23° C (from the start of mixing):
about 2 minutes

3. Pulp capping
Place rubber dam and prepare the cavity. Rinse the cavity and exposed pulpal areas with a suitable disinfectant. Mix a capsule of Harvard MTA-CAP OptiCaps® as described under item 1 and extrude it on a glass plate. With a suitable instrument apply a small amount of Harvard MTA-CAP over the exposed pulp and remove excess moisture with a cotton pellet. At the earliest **3 minutes** after application of Harvard MTA-CAP place a small amount of a flowable light cure liner (e.g. Harvard IonoLine) and light cure. Etch the remaining cavity walls according to the total-etch-technique with Harvard Etch and apply a suitable bonding agent (e.g. Harvard Bond TE Bond) according to the corresponding instructions. Place a light cure composite (e.g. Harvard PremiumFill) according to the instructions and light cure. Check pulp vitality and status regularly.

Additional Notes/Warnings
• In the first hour after application handle the placed MTA cement carefully.
• Store Harvard MTA-CAP OptiCaps® in the sealed packaging in a dry place prior to use.
• Harvard MTA-CAP can cause discoloration.
• Keep away from children!

Composition
MTA, radiopaque mineral oxide

Storage
Store Harvard MTA-CAP at 2 – 25 °C (36 °F – 77 °F) in the original seal-pack and protected from moisture. Do not use after expiry date. Harvard MTA-CAP OptiCaps® are for single use only.

Disposal
Disposal of the product according to local authority regulations.

Reporting obligation
Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note
The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Warranty
Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbHs sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability
Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order.-No.	Article
7081503	Harvard MTA-CAP 2 OptiCaps® ea. 0.25 g, ea. packed in an aluminium pouch
7092000	Harvard Applier OptiCaps®

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