Harvard

MTA-Repair OptiCaps®

(MTA XR)

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Instructions for Use **EN**
Bioactive MTA cement for filling of root perforations, extra radiopacity

Product description
Harvard MTA-Repair OptiCaps® is a biocompatible endodontic cement in OptiCaps®. Due to it's firm consistency in combination with a comfortable setting time Harvard MTA-Repair is the ideal choice for **repair of root perforations**. The powder consists of very fine hydrophilic particles of several hydraulic mineral oxides. After contact with the liquid it forms a gel that hardens to an impermeable barrier. Harvard MTA-Repair OptiCaps® are easily activated and the content of the capsule is easily ejected with the Harvard Applier OptiCaps®. Mixing of Harvard MTA-Repair OptiCaps® (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/minute.

Indications/Intended use

- Repair of root perforations

Performance features
The performance features of the product meet the requirements of the intended use.

Contraindications
None known

Patient target group
Persons who are treated during a dental procedure.

Intended users
This medical device should only be used by a

professionally trained dental practitioner.

Application

1. Activation and mixing
See "Instructions for the activation and mixing of Harvard OptiCaps®"
Mixing time for the Harvard MTA-Repair OptiCaps® is 30 seconds.
Attention:
Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting and lag times may impair or prevent application of the material. The material must be extruded within **10 seconds** after the end of mixing. To prevent dehydration during setting, apply Harvard MTA-Repair immediately after mixing.

2. Working time
Working time at 23° C (from the start of mixing):
about **2 minutes**

3. Repair of root perforations
Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points and isolate the perforation. Fill the root canal apical of the perforation completely with a suitable root canal filling material. Mix a capsule of Harvard MTA-Repair OptiCaps® as described under item 1 and extrude it on a glass plate. Apply Harvard MTA-Repair with suitable instruments into the perforation site and condense it. Check the position of Harvard MTA-Repair in the root canal by an X-ray. If an adequate barrier has not been created, rinse the MTA cement out of the canal and repeat the procedure. Remove excess moisture with a cotton pellet or a paper point. Place a cotton pellet in the access to the root canal and apply a temporary filling material. Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling. Both options can be done at the earliest **5 minutes** after placement of the Harvard MTA-Repair. Harvard MTA-Repair repair material remains as a permanent part of the root canal filling.

Additional Notes/Warnings

- In the first hour after application handle the placed MTA cement carefully.
- Store Harvard MTA-Repair OptiCaps® in the sealed packaging in a dry place prior to use.
- Harvard MTA-Repair can cause discoloration.
- Keep away from children!

Composition
MTA, radiopaque mineral oxide

Storage
Store Harvard MTA-Repair at 2 – 25 °C (36 °F – 77 °F) in the original seal-pack and protected from moisture. Do not use after expiry date. Harvard MTA-Repair OptiCaps® are for single use only.

Disposal
Disposal of the product according to local authority regulations.

Reporting obligation
Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note
The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (<https://ec.europa.eu/tools/eudamed>).

Warranty
Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbHs sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability
Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order.-No.	Article
7081505	Harvard MTA-Repair 2 OptiCaps® ea. 0.25 g, ea. packed in an aluminium pouch
7092000	Harvard Applier OptiCaps®

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