# Harvard MTA Universal OptiCaps®

## Harvard MTA Universal OptiCaps®

### Gebruiksaanweisung DE

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## Instructions for Use EN

### Endodontic repair cement based on mineral oxides

#### Product description

**Harvard MTA Universal OptiCaps®** is a biocompatible, endodontic repair cement in capsules. The powder consists of fine hydrophilic particles of various mineral oxides.

After contact with the liquid it forms a gel that hardens to an impermeable barrier.

Harvard MTA Universal OptiCaps® are easily activated and the content of the capsule is easily ejected with the Harvard Applier OptiCaps®.

Mixing of Harvard MTA Universal OptiCaps® (mixing time 30 seconds) is achieved by a high frequency mixer with about 4,300 oscillations/minute.

#### Indications/Intended use

- Repair of root perforations
- Root-end filling (retrograde)
- Pulp capping (direct and indirect)
- Apexification / Root-end filling (orthograde)
- Repair of internal resorption
- Pulpotomy

#### Performance features

The performance features of the product meet the requirements of the intended use.

#### Contraindications

None known.

#### Patient target group

Persons who are treated during a dental procedure.

#### Intended users

This medical device should only be used by a professionally trained dental practitioner.

#### Application

##### 1. Activation and mixing

See "Instructions for the activation and mixing of Harvard MTA Universal OptiCaps®"

**Mixing time for the Harvard MTA Universal OptiCaps® is 30 seconds.**

#### Attention:

Avoid lag times between the processes of activation, mixing and application as the material is in the process of setting and lag times may impair or prevent application of the material.

The material must be extruded within **10 seconds** after the end of mixing.

To prevent dehydration during setting, apply Harvard MTA immediately after mixing.

#### 2. Working time

Working time at 23° C (from the start of mixing): **about 2 minutes**

#### 3. Applications

##### 3.1. Repair of root perforations

Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points and isolate the perforation.

Fill the root canal apical of the perforation completely with a suitable root canal filling material.

Mix a capsule of Harvard MTA Universal OptiCaps® as described under item 1 and extrude it on a glass plate.

Apply Harvard MTA with suitable instruments into the perforation site and condense it.

Check the position of Harvard MTA in the root canal by an X-ray. If an adequate barrier has not been created, rinse the MTA cement out of the canal and repeat the procedure.

Remove excess moisture with a cotton pellet or a paper point. Place a moist cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

Both options can be done at the earliest **5 minutes** after placement of the Harvard MTA.

Harvard MTA repair material remains as a permanent part of the root canal filling.

##### 3.2. Root-end filling (retrograde)

Create a surgical access to the root-end and resect the root. Prepare an apical cavity to a depth of 3-5 mm.

Isolate the area and dry the root-end cavity with paper points. Achieve hemostasis with suitable methods.

Mix a capsule of Harvard MTA Universal OptiCaps® as described under item 1 and extrude it on a glass plate.

Apply Harvard MTA with suitable instruments and condense it using a small pluggler.

Remove excess cement and clean the surface of the root with a moist piece of gauze.

Confirm placement of Harvard MTA repair material with an X-ray. Harvard MTA remains as a permanent part of the root canal filling.

##### 3.3. Pulp capping

Place rubber dam and prepare the cavity. Rinse the cavity and exposed pulpal areas with a suitable disinfectant.

Mix a capsule of Harvard MTA Universal OptiCaps® as described under item 1 and extrude it on a glass plate.

With a suitable instrument apply a small amount of Harvard MTA over the exposed pulp and remove excess moisture with a cotton pellet.

At the earliest **5 minutes** after application of Harvard MTA place a small amount of a flowable light cure liner (e.g. Harvard IonoLine) and light cure.

Etch the remaining cavity walls according to the total-etch-technique with Harvard Etch and apply a suitable bonding agent (e.g. Harvard Bond TE Bond) according to the corresponding instructions.

Place a light cure composite (e.g. Harvard PremiumFill) according to the instructions and light cure.

Check pulp vitality and status regularly.

#### 3.4. Apexification / Root-end filling (orthograde)

Place rubber dam and clean the root canal system using intracanal instruments and irrigate with NaOCl. Dry the root canal with paper points.

For disinfection place calcium hydroxide paste in the root canal and seal the access opening with a temporary filling material.

After one week remove the calcium hydroxide paste from the root canal system.

Mix a capsule of Harvard MTA Universal OptiCaps® as described under item 1 and extrude it on a glass plate.

With a suitable instrument apply a small amount of Harvard MTA into the apical region and condense it. Create a 3 – 5 mm barrier of MTA cement.

Check the position of Harvard MTA by an X-ray. If an adequate barrier has not been created, rinse the MTA cement out of the canal and repeat the procedure.

At the earliest **5 minutes** after application of the Harvard MTA place a cotton pellet in the access to the root canal and apply a temporary filling material.

Alternatively seal the root canal with a suitable root canal filling material and seal the cavity with a tight filling.

Harvard MTA remains a permanent part of the root canal filling.

#### Additional informations/Warnings

- In the first hour after application handle the placed MTA cement carefully.
- Store Harvard MTA Universal OptiCaps® in the sealed packaging in a dry place prior to use.
- Harvard MTA can cause discoloration.
- Keep away from children!

#### Composition

MTA, bismuth oxide

#### Storage

Store Harvard MTA at 10-25 °C (50-77 °F) in the original seal-pack and protected from moisture. Do not use after expiry date.

Harvard MTA Universal OptiCaps® are for single use only.

#### Disposal

Disposal of the product according to local authority regulations.

#### Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

#### Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>)

#### Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

#### Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order.-No.	Article
7081502	Harvard MTA Universal OptiCaps® 2 OptiCaps® ea. 0.3 g, ea. packed in an aluminium pouch
7092000	Harvard Applier OptiCaps®

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