Harvard SuperSnap® Light Quick Snap

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Instructions for Use EN

Vinylpolysiloxane impression material, light bodied

Product description

Harvard SuperSnap® Light Quick Snap is a super hydrophilic low viscosity impression material. It is characterized by a long working time in conjunction with a short intraoral setting time (super snapset).

Harvard SuperSnap® Light is a precision impression material based on addition curing vinyl silicones delivered in automix-cartridges. It can be extruded directly onto the teeth or the tray.

Harvard SuperSnap® Light is thixotropic. That means it stays in position in the tray, but under pressure it flows. The product shows strong hydrophilic properties. Therefore it improves the wettability and yields an exact reproduction of subgingival areas and other moist surfaces.

Harvard SuperSnap® Light meets the requirements of ISO 4823, type 3 (light-bodied consistency).

Indications/Intended use

- Correction impression technique (2 step/2 phase) with Harvard SuperSnap® Putty
- Double-mix technique (1 step/2 phase) with Harvard SuperSnap® Putty Soft Quick Snap, Harvard SuperSnap® Putty Soft MaxiMix Quick Snap or Harvard SuperSnap® Heavy Quick Snap

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

None known

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Application

1. Preparation of the automix-cartridge

Harvard Dental International GmbH only recommends for Harvard SuperSnap® Light the use of mixing tips type Harvard Auto-T 1:1 O-yellow.

Mount the cartridge into the application gun. Remove and dispose the closure-cap. For exact flow control extrude slightly material until uniform flow from both orifices is achieved. Install the mixing tip onto the cartridge and turn it 90° clockwise until it locks in position and extrude needed material. For direct extruding around the teeth use Intra Tips yellow.

Note:

Store used cartridge with fixed used mixing tip. **Do not re-use the cartridge cap.** Closing the cartridge by turning on the cap between use could cause cross contamination of components and thus premature setting. Leave the used mixing tip on the cartridge. It serves as a cap.

2. Impression taking

Two-step (correction) impression technique

Cut drain channels into the pre-impression. Load the auto mixed Harvard SuperSnap® Light Quick Snap onto the tray with the preliminary impression and/or directly around the teeth. Insert the tray not later than **1:10 minutes** from start of mixing (working time at 23°C or 74°F) into the mouth. Maintain position for at least **1:15 minutes**.

Double-mix technique

Load the impression tray with the tray impression material (with Harvard SuperSnap® Putty Soft Quick Snap, Harvard SuperSnap® Putty Soft MaxiMix Quick Snap or Harvard SuperSnap® Heavy Quick Snap) and apply Harvard SuperSnap® Light Quick Snap around the preparation. Tray loading and syringing around the preparation have to be completed at the same time.

Insert the tray not later than **1:10 minutes** from start of mixing (working time at 23°C or 74°F) into the mouth. Maintain position for at least **1:15 minutes**.

3. Casting

The impression should be stored under normal room conditions (max. 25°C or 77°F). The impression can be poured 30 minutes after removal from the mouth. No other time limitation does exist.

4. Electroplating

All impressions can be copper- or silverplated.

5. Disinfection

Impressions can be disinfected with disinfection solutions specific for impression materials based on addition curing silicones. Refer to the manufacturer instructions.

Storage

Do not store above 25 °C (77 °F). Protect from direct sunlight. Do not use after expiry date.

Additional notes:

- Do not use with standard (condensation reaction type) silicone impression materials.
- Do not use Latex or sulfur-containing polymeric gloves, they may interfere with the setting reaction of the impression material.
- When using the double-mix technique (1 step/2 phase), make sure that both materials have the same setting time in mouth. The working time is determined by the material with the shorter working time.
- Keep away from children!

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7083801	Harvard SuperSnap® Light, Quick Snap, orange, 2 automix cartridges ea. 50 ml, 12 mixing tips 1:1, 4 Intra Tips
7098000	Harvard Auto-T 1:1 O-yellow, refill bag with 50 mixing tips for Light/Medium – Regular/Quick
7083619	Harvard Intra Tips yellow
7095100	refill bag with 50 intra tips. Harvard Dispenser Automix 1:1

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