Harvard TEMP C&B Pro

V03_MDR
2024-10

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Instructions for Use EN

Temporary crown and bridge material, 4:1

Product description

Harvard TEMP C&B Pro is a self cure two component composite in cartridges based on multifunctional methacrylic esters. It is used for creating short- and long-term temporary crowns or bridges, inlays, onlays and veneers. Because of its flexibility, the material is also suitable for longer bridges. Harvard TEMP C&B Pro shows fluorescence in UV-light.

Harvard TEMP C&B Pro is free of methyl methacrylate. Its temperature derived while curing is lower than 40 °C / 104 °F. As a temporary crown or bridge it protects the prepared teeth against external influences and preserves the occlusion. It shows high flexural strength, good abrasion resistance and low polymerization shrinkage. It fits perfectly. Showing good polishability, good color stability and fluorescence its a high esthetic material.

Defects on the temporary made with Harvard TEMP C&B Pro can be easily corrected, as already hardened material aggregates freshly extruded material and light curing composites.

Indications/Intended use

- Fabrication of temporary crowns, bridges, inlays, onlays and veneers

Performance features

The performance features of the product meet the requirements of the intended use.

Contraindications

Do not use for patients with allergic reactions against acrylates. Harvard TEMP C&B Pro should not be used any more, if allergic reactions are observed.

Patient target group

Persons who are treated during a dental procedure.

Intended users

This medical device should only be used by a professionally trained dental practitioner.

Application

1. Impression taking

Before preparation or extraction of the tooth or teeth make a situation impression with addition curing silicones (long storage stability) or with alginates. For better stability of the temporary carve out interdental areas. In molar areas with teeth absent it may be necessary to cut a groove in the impression between the abutments to create a bridge-like connection between the tooth units.

Note:

Block out undercuts and if necessary cut grooves into the impression.

2. Preparing of the Automix-cartridge

Remove the cap of the Automix-cartridge and throw it away (**do not use it again!**). Attach a 4:1 mixing tip. Make sure that the guidance of the Automix-cartridge is aligned with that of the mixing tip and turn the tip 90° clockwise until it locks in position. When the application gun is loaded with the prepared cartridge it is ready for application.

Note:

Discard the initial extrusion from the mixing tip (about the size of a pea). Then the following mix will be perfect. This must be done for each new mix.

Leave the used mixing tip on the cartridge. It serves as a cap.

3. Forming of the temporary crowns / bridges

Harvard TEMP C&B Pro mixes automatically when dispensed with slight and even pressure directly into the situation impression made before. Filling should occur from bottom upward to prevent voids. **Working time of Harvard TEMP C&B Pro is 50 seconds (at 23 °C/74 °F).**

3.1. Application in the mouth

Load the situation impression with Harvard TEMP C&B Pro. Seat the impression within the working time (**50 seconds** at 23 °C/74 °F) onto the prepared teeth. After **1-2 minutes** (setting time in mouth at 37 °C/98 °F) the material shows a hardened but still elastic condition and can be removed from the teeth together with the situation impression.

3.2. Application on the model

Load the situation impression with Harvard TEMP C&B Pro. Seat the impression within the working time (**50 seconds** at 23 °C/74 °F) onto the prepared areas of the model. After **3-4 minutes** (setting time at 23 °C/74 °F) the material shows a hardened but still elastic condition and can be removed from the model together with the situation impression.

Note:

The setting reaction has to be checked with excess material intra-orally (e.g. with a scaler) resp. on the model. The oral temperature has a significant effect on the setting reaction and the temporary can only be removed without destruction during the elastic state.

4. Post-curing and finishing

If possible, leave the temporary in the situation impression during post-curing. Optimally the temporary is post-cured in warm water (45-55 °C/ 113-131 °F, e.g. in a hot cure polymerization device) until reaching its final hardness (approx. **4:30** minutes from start of mixing). Post-curing at room temperature is completed after about **6 minutes** from start of mixing. When the cure is completed remove the temporary. Afterwards remove excess material and proximal undercuts. Then the temporary can be refined with rotary instruments and polished to high gloss.

Do not breathe polishing dust; use suitable mouth protective device, safety glasses and aspiration!

Note:

The oxygen inhibited smear layer on the surface caused by the acrylic system can not be fully avoided and should be removed. It can easily be removed by alcohol or other suitable solvents.

5. Cementing of the temporary

Harvard TEMP C&B Pro temporary crowns or bridges preferably should be cemented preferably with an eugenol free temporary cement (e.g. Harvard TEMP Cem). In case of using eugenol containing cements it cannot be excluded that later used acrylate based cements could be hindered in setting (inhibition of curing reaction by traces of remained eugenol).

6. Repairs

Harvard TEMP C&B Pro temporaries show high mechanical strength. However, if a temporary breaks the following procedures are recommended:

6.1. Fracture shortly after production

Both ends of the fracture are refixed with freshly extruded Harvard TEMP C&B Pro.

6.2. Fracture of longer existing worn temporaries

Clean and roughen the areas of fracture and provide with some mechanical retention. Thus prepared, join the fractured areas with freshly extruded Harvard TEMP C&B Pro. To enhance the curing place for some minutes in warm water (50 °C / 122 °F).

7. Remarks

Non cured Harvard TEMP C&B Pro can easily be removed with alcohol or other suitable solvents.

8. Process overview

Up to 50 s	impression loading at 23 °C / 73 °F and application in mouth
50 s – 3 min	setting in mouth at 37 °C / 99 °F
3 – 4,5 min	polymerization at ca. 50 °C / 122 °F
From 4,5 min	trimming / polishing at 23 °C / 73 °F

Additional informations/Warnings

- Harvard TEMP C&B Pro is free of methylnmethacrylate but contains other methacrylates.
- Avoid contact with skin, mucous membrane and eyes
- If the material comes into contact with skin, immediately wash with water and soap. If the material comes into contact with eyes, immediately rinse with copious amounts of water and seek medical advice if required.
- For sensitive persons, sensitization to the material cannot be excluded. The material should not be used any more, if allergic reactions are observed.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- Keep away from children!

Composition

Unsaturated esters, multifunctional methacrylates, multifunctional acrylates, malonyl urea derivatives, dental glass, silicon dioxide, catalysts

Disposal

Disposal of the product according to local authority regulations.

Reporting obligation

Serious incidents according to the EU Medical Devices Regulation that have occurred in connection with this medical device must be reported to the manufacturer and the competent authority.

Note

The summary of safety and clinical performance of the medical device can be found in the European database on medical devices (EUDAMED – <https://ec.europa.eu/tools/eudamed>).

Storage

Do not store above 25 °C (77 °F). Do not use after expiration date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

Keep the instructions for use for the duration of the application.

Order No.	Article
7081641	50 ml automix cartridge, ratio 4:1, 10 mixing tips
7081642	Shade A1
7081643	Shade A2
7094000	Shade A3
	Harvard Auto 4:1 / 10:1 S-Blue, refill bag with 50 mixing tips
7095000	Harvard Dispenser Automix 4:1 / 10:1

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