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**Thixotropic gel
containing 37% phosphoric acid**

Harvard Etch is a 37% phosphoric acid gel with an excellent thixotropy for use in the total etch technique.

Indications

Enamel etch or total etch technique for:

- Composite restorations
- Sealing of fissures
- Adhesive cementation of inlays, onlays, crowns and bridges (with suitable composite-cement systems)

Contraindications

The application of Harvard Etch is contraindicated if the recommended application technique is not possible.

Side effects

Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection the cavity floor in cases of deep excavations should be covered with a thin layer of calcium hydroxide material.

Application**1. Preparation of Harvard Etch syringe**

Screw off the cap and substitute it by the application needle. **Store the screw cap!**

Note:

The application needle might be bend in a slight arc, if desired, to facilitate application.

Handling precaution:

Squeeze out the first small drop on a pad to ensure proper handling and prevent "shooting" out of the etching gel by entrapped air.

2. Application of Harvard Etch

Recommended is the total etch technique: Apply Harvard Etch onto the enamel and dentine surfaces beginning with the enamel bevels. **Condition the enamel for at least 15 seconds and the dentine for 15 seconds.** Deciduous teeth are etched correspondingly longer. Rinse for 20 seconds with water. Dry it in a water and oil-free airstream, but **do not** desiccate. A slightly wet dentine surface is important for the function of the bonding agent (e.g. Harvard Bond TE Mono). The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

After use replace disposable application cannula by the screw cap. Discard used application cannula.

3. Application of the bonding agent

Apply a bonding agent (e.g. Harvard Bond TE Mono) with a brush on the etched saliva free enamel and dentine surfaces according to the corresponding instructions.

Warnings

- Harvard Etch contains phosphoric acid and is corrosive. Avoid contact with eyes, mucous membrane and skin (safety glasses are recommended for both patients and dentists). After contact with eyes immediately rinse with copious amounts of water for at least 15 minutes and consult a physician/ophthalmologist. After skin contact immediately rinse with soap and water.

- Avoid contact with gingiva, certain liners (based on glass ionomers) and adjacent teeth.

Storage information

Storage temperature 4-25°C (39-77°F). Do not use after expiry date.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbHs sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!**Store product out of reach of children!**

Keep the instructions for use for the duration of the application.

Order-No.	Article
	Harvard Etch
7054000	2ml syringe, 3 needle tips
7095126	Harvard NeedleTips H25
	refill bag with 50 needle tips

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