Harvard Restore System Kit Plus:

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7083217

Harvard Restore System Kit Plus:

- **Harvard Restore** is a light cure microhybrid composite for anterior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05-1.5 µm. The total filler content is 81 % by weight and 65 % by volume.
- **Harvard RestoreFlow** is a visible light cure easy flowable microhybrid composite for anterior and limited posterior restorations and for sealing. Harvard RestoreFlow is based on methacrylate-resin and inorganic filler particles of 0.05-1 µm. The filler content is 61 % by weight and 41 % by volume.
- **Harvard Self-Bond** is a simple to use self-etching light curing one component adhesive. It ensures a strong bonding of composites to enamel and dentin. Harvard Self-Bond is compatible with all current brands of visible light cure methacrylate based composite restorative materials.

Indications

Suitable for all cavity classes.

Contraindications

The placement of Harvard Restore restorations is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard Restore if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Harvard Restore may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite come into contact with materials containing eugenol.

Application

1. Preparation of tooth and color matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching. Ascertain the tooth shade while teeth are still moist and select the appropriate Harvard Restore shade. Harvard Restore is shaded according to VITA®-Shades.

2. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard CalciumHydroxide).

3. Application of Harvard Self-Bond

Apply Harvard Self-Bond with a brush onto the **slightly moist enamel and dentin surfaces** for at least 15 seconds with agitation.

The material should build a homogeneous layer. Air thin gently for 5-10 seconds to remove the volatile components and to disperse the adhesive. Take care, that the adhesive is not removed from the tooth surface by a strong air stream.

Do not touch the treated surface!

Then light cure with a suitable polymerization unit (wavelength range 400–500 nm, light intensity of at least 500mW/cm²) for 20 seconds.

4. Application and curing of Harvard Restore

Apply Harvard Restore in the selected shade. Place it by using suitable instruments. Harvard Restore is easy to model. Matrix strips may be used.

In case of deep and extended cavities best results are obtained with a built-up of a base liner using a light cure resin modified cement (e.g. Harvard Ionoglas Fill LC) or a light cure flowable composite (e.g. Harvard RestoreFlow), followed by the application of Harvard Restore. In case of small cavities Harvard Restore can be applied directly.

For an optimum result apply Harvard Restore in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit (wavelength range 400 – 500 nm) with a light intensity of at least 1000mW/cm², cure each increment as follows:

Lighter shades (e.g. A1, A2, C2) **20 seconds**
Darker shades (e.g. A3.5, B3) **30 seconds**
Hold the light emission window as close as possible to the filling material.

5. Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burs. Polish to high gloss with polishing discs or polishing brushes.

6. Additional informations

- Do not use any resin to adjust viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard Restore and Harvard Self-Bond may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the product for user's application. If this product is defective within the warranty period, your exclusively remedy and Harvard Dental International's sole obligation shall be repair or replacement of the Harvard Dental International product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Harvard Restore conforms to the ISO 4049.

Keep the instructions for use for the duration of the application.

Order No. Article

- System Kit Plus:
4 syringe Harvard Restore ea. 4g
3 ml syringe Harvard RestoreFlow A2
5 ml bottle Harvard Self-Bond
5 needle tips, 50 x Microbrush[®],
5 mixing wells
Harvard Restore
Shades A1, A2, A3, A3.5

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