Harvard Restore System Kit TE Mono

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Harvard Restore System Kit TE Mono:

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Instructions for Use EN

Harvard Restore System Kit TE Mono:

- **Harvard Restore** is a light cure microhybrid composite for anterior and posterior restorations. It is based on BIS-GMA-resin and inorganic filler particles of 0.05-1.5 µm. The total filler content is 81 % by weight and 65 % by volume.
- **Harvard Etch** is a 37% phosphoric acid gel with an excellent thixotropy for use in the total etch technique.
- **Harvard Bond TE Mono** is a simple to use light cure one component bonding agent for the total etch technique with a phosphoric acid etching gel ("Etch and Rinse"). It is designed for strong bonding of light cure composites to etched enamel and dentine and to nonprecious and precious metals. Good, long lasting adhesive strength and good biocompatibility are attained by polycarboxylic acids. Harvard Bond TE Mono is compatible with all current brands of visible light cure methacrylate-based composite materials. It is ethanol-based and hydrophilic. Harvard Bond TE Mono can be used on slightly moist dentin surfaces (wet bonding technique).

Indications

Suitable for all cavity classes.

Contraindications

The placement of Harvard Restore restorations is contraindicated if a dry working area or the recommended application technique are not possible. Also do not use Harvard Restore if the patient is known to be allergic to any of the ingredients.

Side effects

In singular cases, Harvard Restore may cause a sensitizing reaction in patients with a hypersensitivity to any of the ingredients. In these cases, the material should not be used. Irritations resulting from direct contact with the pulp cannot be ruled out. Therefore for pulp protection areas close to the pulp should be covered with a thin layer of calcium hydroxide material.

Incompatibility with other materials

Do not use in combination with substances containing eugenol because eugenol inhibits the polymerization of the composite. Neither store the composite material in proximity of eugenol containing products, nor let the composite allow coming into contact with materials containing eugenol.

Application

1. Preparation of tooth and color matching

Clean the tooth with a fluoride-free polishing paste (e.g. flour of pumice with water) prior to preparation and color matching.

Ascertain the tooth shade while teeth are still moist and select the appropriate Harvard Restore shade.

Harvard Restore is shaded according to VITA®-Shades.

2. Cavity preparation

After isolation (best with a rubber dam) prepare the cavity with minimal tooth reduction. For pulp protection areas close to the pulp should be covered with a thin layer of calciumhydroxide material (e.g. Harvard Calcium-Hydroxide).

3. Application of Harvard Etch

Recommended is the total etch technique: Apply Harvard Etch onto the enamel and dentine surfaces beginning with

the enamel bevels. Condition the enamel for at least **15 seconds** and the dentine for **15 seconds** or less. Deciduous teeth are etched correspondingly longer. Rinse for **20 seconds** with water. Dry it in a water and oil-free airstream, but do not desiccate. A slightly wet dentine surface is important for the function of the bonding agent (e.g. Harvard Bond TE Mono). The etched enamel bevel should have a chalky white appearance.

Etching precaution:

It is essential, that etched areas are not contaminated by anything. If contamination occurs, re-etch, rinse with water and dry as above. Avoid etching gel contact with oral soft tissues, eyes and skin. If accidental contact occurs, flush immediately with copious amounts of water.

After use replace disposable application cannula by the screw cap. Discard used application cannula.

4. Application of Harvard Bond TE Mono

Apply Harvard Bond TE Mono generously with a brush onto the enamel and dentin surfaces for **30 seconds** with agitation. The material should build a homogeneous layer. Remove excess material carefully. Dry cautiously with oil free air for about **15 seconds** to remove all volatile components and to disperse the adhesive to an even layer. Do not desiccate the dentine.

Then light cure with a suitable polymerization unit (wavelength range 400–500 nm, light intensity of at least 500mW/cm²) for **20 seconds** before application of a second layer of Harvard Bond TE Mono.

Apply again the Harvard Bond TE Mono generously with a brush onto the adhesive surfaces as described above before placement of a light cure composite.

Do not touch the treated surface!

5. Application and curing of Harvard Restore

Apply Harvard Restore in the selected shade. Place it by using suitable instruments. Harvard Restore is easy to model. Transparent matrix strips may be used. In case of deep and extended cavities best results are obtained with a build-up of a base liner using a light cure resin modified cement (e.g. Harvard Ionoglas Fill LC) or a light cure flowable composite (e.g. Harvard RestoreFlow), followed by the application of Harvard Restore. In case of small cavities Harvard Restore can be applied directly.

For an optimum result apply Harvard Restore in **layers of max. 2 mm**. Light cure each increment separately.

By using a polymerization unit (wavelength range 400 – 500 nm) with a **light intensity of at least 1000mW/cm²**, cure each increment as follows:
Lighter shades (e.g. A1, A2, C2) **20 seconds**
Dark shades (e.g. A3.5, B3) **30 seconds**
Hold the light emission window as close as possible to the filling material.

6. Finishing

Depending on the shape of the area contour with finishing diamonds, flexible grinding discs or carbide burrs. Polish to high gloss with polishing discs or polishing brushes.

7. Warnings

- Harvard Etch contains phosphoric acid and is corrosive. Avoid contact with eyes, mucous membrane and skin (safety glasses are recommended for both patients and dentists). After contact with eyes immediately rinse with copious amounts of water for at least 15 minutes and consult a physician/ophthalmologist. After skin contact immediately rinse with soap and water.
- Avoid contact with gingiva, certain liners (based on glass ionomers) and adjacent teeth.

8. Additional informations

- Do not use any resin to adjust the viscosity of composite restorative material.
- Avoid contact with skin, mucous membrane and eyes.
- Unpolymerized Harvard Restore and Harvard Bond TE Mono may have an irritant effect and can lead to sensitization against methacrylates.
- Commercial medical gloves do not protect against the sensitizing effect of methacrylates.
- VITA® is a registered trade mark of the VITA-Zahnfabrik, Bad Säckingen, Germany.

Storage information

Storage temperature 4-25°C (39-77°F). Protect from direct sunlight. Do not use after expiry date. Close syringe immediately after use to avoid exposure to light.

Warranty

Harvard Dental International GmbH warrants this product will be free from defects in material and manufacture. Harvard Dental International GmbH makes no other warranties including any implied warranty of merchantability or fitness for a particular purpose. User is responsible for determining the suitability of the

product for user's application. If this product is defective within the warranty period, your exclusive remedy and Harvard Dental International GmbH's sole obligation shall be repair or replacement of the Harvard Dental International GmbH product.

Limitation of Liability

Except where prohibited by law, Harvard Dental International GmbH will not be liable for any loss or damage arising from this product, whether direct, indirect, special, incidental or consequential, regardless of the theory based asserted, including warranty, contract, negligence or strict liability.

For dental use only!

Store product out of reach of children!

Harvard Restore conforms to the ISO 4049.

Keep the instructions for use for the duration of the application.

Order-No.	Article
7083207	Harvard Restore System Kit TE Mono 4 syringe ea. 4g shade A1, A2, A3, A3.5 5 ml bottle Harvard Bond TE Mono 3 ml syringe Harvard Etch 5 needle tips, 50 x Microbrush®, 5 mixing wells

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