

		suitable for all dental conditioning techniques: with or without etching gel, wet or dry (1 component, 1 step)							
		Activator for safe adhesion for all applications where the light-curing of Harvard InterLock® ONE can not be guaranteed							
		Dual-curing, self-etching adhesive (2 components, 1 step)							
		Self-etching, light-curing adhesive (1 component, 1 step)							
		Very strong, self-etching, light-curing adhesive (2 components, 2 steps)							
		Self-etching, light-curing adhesive (1 component, 1 step)							
		Light-curing adhesive for the total-etch technique with one thixotropic etching gel							
Dual-Curing		X	X						
Light-Curing		X		X	X	X	X		X
Self Etch		X	X	X	X	X			
Total-Etch-Technique		X							X
Universal Adhesive	X								
	Harvard InterLock® ONE								
	Harvard InterLock® ONE Self Cure Activator*								
	Harvard Bond SE Dual								
	Harvard Self-Bond								
	Harvard InterLock®								
	Harvard Bond SE Mono								
	Harvard Bond TE Mono + Harvard Etch								

* in combination with Harvard InterLock® ONE

The Harvard Bonding Competence

the overview:

Harvard InterLock® ONE Universal Adhesive



Harvard InterLock® ONE Self Cure Activator



Harvard Bond SE Dual



Harvard Self-Bond



Harvard InterLock®



Harvard Bond SE Mono



Harvard Bond TE Mono



Harvard Etch



Really universal: suitable for all

teeth conditioning techniques:

as "Self Etch" adhesive, as "Total Etch" adhesive (which etching step), wet or dry

- Only one thin layer necessary
- High & reproducible bond strength to enamel and dentine even under long term stress
- Effective marginal seal for low risk of postoperative sensitivity
- Excellent performance also with self or dual cure composites when light curing is possible, if not see e.g. for adhesive post cementing in the root canal mix additional with **Harvard InterLock ONE Self Cure Activator****
- Convenient application by bottle (precise dropper) or SingleDose



Self Cure Activator

**

- Activator for safe adhesion in indications, where light curing of Harvard InterLock® ONE cannot be ensured (e.g. in the root canal)
- Ideal for cementing of posts with dual and self cure composites

Article information

	Order no.
Harvard InterLock® ONE Universal Adhesive 5 ml bottle	7083611
Microbrush® refill bag with 50 applicators	7095156
Harvard InterLock® ONE SingleDose Universal Adhesive 50 x 0.1 ml SingleDose	7083613
Harvard InterLock® ONE Self Cure Activator 5 ml bottle	7083612
Harvard Bond SE Dual 2 x 5 ml bottle, 50 Microbrush®, 5 mixing wells	7083602
2 x 5 ml bottle (Refill)	7083605
Microbrush® refill bag with 50 applicators	7095156
Harvard Self-Bond 5 ml bottle	7083218
Harvard InterLock® 2 x 5 ml bottle, 50 Microbrush®, 5 mixing wells	7083603
Microbrush® refill bag with 50 applicators	7095156
Harvard Bond SE Mono Kit 5 ml bottle, 25 Microbrush®, 5 mixing wells	7083604
Harvard Bond SE Mono Refill 5 ml bottle	7083601
Harvard Bond SE Mono SingleDose Kit 50 SingleDose ea. 0.15 ml, 50 Microbrush®	7083609
Microbrush® refill bag with 50 applicators	7095156
Harvard Bond TE Mono 5 ml bottle	7083607
Harvard Etch 2 ml syringe, 3 needle tips	7054000
Harvard NeedleTips H25 refill bag with 50 needle tips	7095126



HARVARD®

Marke und Qualität seit 1892

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Phone: + 49 (0) 30/99 28 978-0
Fax: + 49 (0) 30/99 28 978-19
info@harvard-dental-international.de

10996_EN_09/18

HARVARD®

Harvard InterLock® ONE

Universal Adhesive

NEW



5 ml



0.1 ml



Please rate this product.*

EN

1. Light Cure Composites, Bondings and Liner



Harvard UltraFill

Harvard PremiumFill

Harvard UltraFlow



Harvard InterLock® ONE Universal Adhesive



Harvard Restore



Harvard RestoreFlow



Harvard InterLock®



Harvard Bond SE Mono



Harvard IonoLine



Harvard Self-Bond



Harvard Bond TE Mono



Harvard Etch



Harvard CalciumHydroxide LC



Harvard ZirconCore



Harvard InterLock® ONE Self Cure Activator



Harvard MTA XR Flow Fast OptiCaps®



Harvard MTA XR Flow EWT OptiCaps® & EndoDirect syringe



Harvard Bond SE Dual



Harvard LuteCem SE



Harvard MTA XR Fast OptiCaps®



Harvard MTA XR OptiCaps®



Harvard MTA (Universal) OptiCaps®



Harvard MTA (Universal) HandMix



Harvard TEMP C&B



Harvard TEMP Glaze LC



Harvard TEMP Cem Esthetic



Harvard Implant Semi-Permanent




Harvard TEMP C&B Pro



Harvard TEMP Cem EasyDose®



Harvard Phosphat Cement



Harvard Cement OptiCaps®



Harvard Polycarboxylat Cement

5.2 Glass Ionomer Cements for Restorations



Harvard Ionocoat LC



Harvard Ionoresin Prime LC



Harvard Ionoglas Fill



Harvard Ionoglas Fill LC



Harvard Ionoresin Fill Extra



Harvard Ionoresin Cem Extra



Harvard Ionoresin Cem



Harvard PremiumSil Putty Soft/Heavy Fast

5.3 Glass Ionomer Cements for Luting



Harvard TransMatrix



Harvard PremiumSil MaxMix Putty Soft/Heavy Fast



Harvard PremiumSil Heavy

6. Bite Registration and Impression Materials



Harvard TraySive



Harvard Bite



Harvard PremiumSil Putty Soft/Heavy Fast



Harvard PremiumSil Putty Soft/Heavy Fast

Harvard PremiumSil/Light/Fast, Medium/Medium Fast, Mono/Mono Fast, Heavy/Heavy Fast



*Great that you have tested our product. Did it meet your satisfaction? Your opinion is very important to us; please use the QR code or our website www.harvard-dental.de/feedback to provide us with your much valued feedback. Product tests are worth it. As a gesture of thanks we will send you another product sample from our range. Thank you for your time.

Please rate this product.*

