

Harvard UltraFill concept of opacities.

Overview:




Pure esthetics in a new shine.

1. Multi Opacity layering technique for specialists

Perfect esthetic results with Harvard UltraFill Multi Opacity. With this opalescent high performance composite it is possible to create a perfect solution for challenging restorations. This concept offers a huge variety of Vita®- and special shades.

(Layers of dentin, enamel, incisal and more translucent shades).

Multi Opacity layering technique for specialists



Easy and efficient.

2. Universal Opacity

Easy layering with Harvard UltraFill Universal. High performance results in less time.

(Easy technique – medium opacity).

2 different concepts of opacity

Universal Opacity



A new high performance composite with excellent properties

- Superior handling properties
- Easy to work with
- Easy polishability, permanent high gloss
- Opalescence, tooth-like fluorescence, cameleon effect
- Very good radiopacity
- Low polymerization shrinkage
- High abrasion resistance
- Optimal physical properties
- Inorganic filler smaller than 0,2 µm
- Available in a new 3 g syringe and in OptiTips®

Article information

	Order no. 3 g syringe	Order no. 20 x 0.25 g OptiTips®
Shade A1 E	7131111	7131211
Shade A1 D	7121111	7121211
Shade A2 E	7131112	7131212
Shade A2 D	7121112	7121212
Shade A3 E	7131113	7131213
Shade A3 D	7121113	7121213
Shade A3.5 E	7131114	7131214
Shade A3.5 D	7121114	7121214
Shade B1 E	7131121	7131221
Shade B1 D	7121121	7121221
Shade B2 E	7131122	7131222
Shade B2 D	7121122	7121222
Shade C2 E	7131132	7131232
Shade C2 D	7121132	7121232
Shade A4 E	7131115	7131215
Shade B3 E	7131123	7131223
Shade C3 E	7131133	7131233
Shade INT D	7121140	7121240
Shade BL E	7131153	7131253
Shade BL 0.5 (E)	7111150	7111250
Shade BL D	7121152	7121252
Shade INC	7141160	7141260
Shade T	7141170	7141270
Shade TX	7141180	7141280

Harvard UltraFill „Multi Opacity“ Core Kit 7186101
14 syringes ea. 3 g, in the shades
A1 E, A2 E, A3 E, A3.5 E, B1 E, B2 E, C2 E
A1 D, A2 D, A3 D, A3.5 D, B1 D, B2 D, C2 D

Harvard UltraFill „Multi Opacity“ Complementary Kit 7186102
10 syringes à 3 g, in the shades
A4 E, B3 E, C3 E, BL 0.5 (E), BL E, T, TX, INT D, BL D, INC

Harvard UltraFill „Multi Opacity“ Trial Kit S 7176104
3 syringes ea. 3 g, in the shades INC, A2 E, A2 D

Harvard UltraFill „Multi Opacity“ Trial Kit C 7176204
3 x 20 OptiTips® ea. 0.25 g, in the shades INC, A2 E, A2 D

	Order no. 3 g syringe	Order no. 20 x 0.25 g OptiTips®
Shade A1 U	7111111	7111211
Shade A2 U	7111112	7111212
Shade A3 U	7111113	7111213
Shade A3.5 U	7111114	7111214
Shade B2 U	7111122	7111222
Shade C2 U	7111132	7111232
Shade D3 U	7111143	7111243
Shade BL U	7111151	7111251

Harvard UltraFill „Universal Opacity“ IntroKit S 7175100
3 syringes ea. 3 g, in the shades A2 U, A3 U, A3.5 U

Harvard UltraFill „Universal Opacity“ IntroKit C 7175200
3 x 20 OptiTips® ea. 0.25 g, in the shades A2 U, A3 U, A3.5 U

Harvard Applier OptiTips® 7095200

HARVARD®
Marke und Qualität seit 1892

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Tel.: + 49 (0) 30/99 28 978-0
Fax: + 49 (0) 30/99 28 978-19
info@harvard-dental-international.de

10967_EN_05/18

HARVARD®

Harvard UltraFill Ultra Gloss, Opalescent, High Performance

Composite for anterior and posterior restorations.

WORLD NOVELTY



Please rate this product.*

EN

1. Light Cure Composites, Bondings and Liner



WORLD NOVELTY

Harvard UltraFill



Harvard PremiumFill



Harvard UltraFill Flow



Harvard Inter-Lock® ONE Universal Adhesive



Harvard Restore



Harvard RestoreFlow



Harvard Inter-Lock®



Harvard Bond SE Mono



Harvard IonoLine



Harvard Bond TE Mono



Harvard Etch



Harvard CalciumHydroxide LC



Harvard ZirconCore

Harvard Inter-Lock® ONE Self Cure Activator



Harvard MTA XR Flow Fast OptiCaps®



Harvard MTA XR Fast OptiCaps®



Harvard MTA (Universal) OptiCaps®



4. Temporary Materials

Harvard TEMP C&B



Harvard TEMP C&B Pro



Harvard Phosphat Cement



Harvard Ionoglas Fill Extra



Harvard Ionoresin Fill Extra



Harvard Ionoglas Cem Extra



Harvard Bite



Harvard PremiumSII Light



Harvard PremiumSII Medium



Harvard PremiumSII Heavy

2. Dual Cure Core Build-Up and Post Cementing



Harvard LuteCem SE



Harvard Bond SE Dual



Harvard MTA XR Fast OptiCaps®



Harvard CalciumHydroxide



Harvard TEMP Cem Esthetic



Harvard TEMP Cem EasyDose®



Harvard Cement OptiCaps®



Harvard Polycarboxylat Cement



Harvard MTA XR Flow EWT OptiCaps® & EndoDirect syringe



Harvard MTA XR OptiCaps®



Harvard Implant Semi-permanent

5.1 Zinc Oxide Luting Cements



Harvard Phosphat Cement



Harvard Ionoglas Fill Extra



Harvard Ionoresin Fill Extra



Harvard Ionoglas Cem Extra



Harvard Bite



Harvard PremiumSII Light



Harvard PremiumSII Medium



Harvard PremiumSII Heavy

5.2 Glass Ionomer Cements for Restorations



Harvard Ionoglas Fill Extra



Harvard Ionoresin Fill Extra



Harvard Ionoglas Cem Extra



Harvard Bite



Harvard PremiumSII Light



Harvard PremiumSII Medium



Harvard PremiumSII Heavy



Harvard Ionoglas Fill Extra



Harvard Ionoresin Fill Extra



Harvard Ionoglas Cem Extra



Harvard PremiumSII Max:Mix Putty Soft/Heavy Fast



Harvard PremiumSII Heavy



Harvard PremiumSII Putty Soft/Heavy Fast

Harvard PremiumSII Light/Fast, Medium/Medium Fast, Mono/Mono Fast, Heavy/Heavy Fast



Harvard TraySive



Please rate this product.*

*Great that you have tested our product. Did it meet your satisfaction? Your opinion is very important to us; please use the QR code or our website www.harvard-dental.de/feedback to provide us with your much valued feedback. Product tests are worth it. As a gesture of thanks we will send you another product sample from our range. Thank you for your time.