

Harvard **NEW** TEMP Cem Esthetic

Translucent eugenol-free zinc oxide cement for temporary luting.

Esthetic translucent cement.

- No shine-through of core
- Easy removal of excess material
- Easy removal of crowns and bridges
- After removal no residual cement remains on the core
- Low film thickness for a perfect fit
- Eugenol-free
- Comfortable handling and time saving thanks minimix syringe


Translucency for high esthetics!



Harvard **NEW** LiquiCord®

Retraction paste for temporary gingiva displacement

- For a dry and expanded sulcus
- For patients: less painfull
- For dentists: time saving, easy to use
- Neutral taste



Optional Accessories:



Harvard **NEW** BioCal®-CAP

Bioactive light cure, resin-modified MTA Cement

The plaster for the pulp.




- Calcium release and a high pH promotes the formation of calcium hydroxyapatite respectively tertiary dentin
- For direct and indirect pulp capping for the treatment of deciduous and permanent teeth
- Creates an environment hostile to bacteria
- Moisture tolerant
- No solubility
- Radiopaque
- Light cure



Harvard **NEW** SuperSnap®

Superhydrophilic and SuperSnap vinylpolysiloxane impression material

- Convenient working time – **shortened** intra oral setting time
- Available in two Snap-Set times: **Regular Snap** und **Quick Snap**
- Comfortable, user- and patient-friendly
- Excellent flow properties
- Intensive colors for excellent legibility



Fast as the fastest!



 **HARVARD®**

 **HARVARD®**

“Find your Harvard Smile.”

Innovative Dental Products 2019. Made in Germany.

New at Harvard

- Harvard SuperSnap®
- Harvard LiquiCord®
- Harvard BioCal®-CAP
- Harvard TEMP Cem Esthetic



Much more than you expect.

Much more than you expect.

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Phone: + 49 30/99 28 978-0, Fax: + 49 30/99 28 978-19,
info@harvard-dental-international.de
www.harvard-dental.de

1892
sell
HARVARD

I. Restorative

1: Composite for Fillings



Harvard UltraFill Flow

Harvard PremiumFill®



Harvard PremiumFlow®

Harvard Restore

Harvard RestoreFlow

2: Bondings / Adhesive



Harvard InterLock® ONE Universal Adhesive

Harvard Bond SE Mono

Harvard Self-Bond



Harvard Bond TE Mono

Harvard Etch



Harvard CalciumHydroxide LC

Harvard IonoLine



Harvard Ionoresin Filli LC

Harvard IonoResin Fill Extra

Harvard IonoResin Fill



Harvard Ionocoat LC

Harvard IonoResin Prime LC



Harvard ZirconCore

Harvard InterLock® ONE Self Cure Activator

Harvard Bond SE Dual

II. Endodontics

6: Pulp Capping



Harvard BioCal®-CAP

7: Root Canal



Harvard MTA XR Flow Fast OptiCaps®

Harvard MTA XR Fast OptiCaps®

Harvard MTA XR OptiCaps®

Harvard MTA XR Flow EWT OptiCaps® & EndoDirect syringe



Harvard MTA (universal) OptiCaps®

Harvard MTA (universal) HandMix

Harvard CalciumHydroxide

III. Prosthetics

8: Impressions



Harvard SuperSnap® Light

Harvard SuperSnap® Medium

Harvard SuperSnap® Mono

Harvard SuperSnap® Heavy



Harvard SuperSnap® MaxiMix Putty Soft / Heavy

Harvard PremiumSil Putty Soft / Heavy Fast

Harvard PremiumSil Putty Soft / Putty Soft Fast

Harvard PremiumSil Heavy / Heavy Fast



Harvard PremiumSil Light / Light Fast

Harvard PremiumSil Medium / Medium Fast

Harvard PremiumSil Mono / Mono Fast

Harvard PremiumSil Heavy / Heavy Fast



Harvard TraySive

Harvard Bite

9: Temporary Materials



Harvard TEMP C&B

Harvard TEMP Glaze LC

Harvard TEMP C&B Pro

Harvard TEMP Cem



Harvard TEMP Cem EasyDose®

Harvard TEMP Cem Esthetic

10: Luting Cements



Harvard Implant Semi-permanent

Harvard LuteCem SE

Harvard Polycarboxylat Cement

Harvard Phosphat Cement

11: Glass Ionomer Cement for Luting



Harvard IonoResin Cem Extra

Harvard IonoResin Cem LC

Harvard IonoResin Cem Extra

Harvard IonoResin Cem

12: Transparent Matrix Material



Harvard TransMatrix



Great that you have tested our product. Did it meet your satisfaction? Your opinion is very important to us. **Product tests are worth it.** As a gesture of thanks we will send you another product sample from our range. **Thank you for your time.**

www.harvard-dental.de/feedback

Please rate this product.