


Harvard PremiumFill ⁺

Overview:



Two different concepts of opacity.

1.) Multi opacity layering technique for specialists

Dentin Opacity


Dentin opacity for multi opacity layering technique.

Enamel Opacity

Enamel opacity for multi opacity layering technique.

Transpa Opacity

Transpa opacity as well as other translucent shades for multi opacity layering technique.



2.) Universal opacity for easy layering technique

Easy technique – **one** medium opacity, suitable for easy cases.



The composite with a certain plus.

Long-lasting, esthetic, superior in handling.

Superior handling properties and non-sticky consistency

- Easy to place and convenient to use
- Highly sculptable
- Excellent high gloss polishability
- Inorganic filler about 0,6 µm

Low polymerization shrinkage thanks to the nano fillers

- Perfect marginal seal

Perfect esthetic restorations

- Tooth-like fluorescence
- Chameleon-like color adaptation
- Large choice of shades with different translucencies

Delivered in 4 g syringes and 0.3 g OptiTips[®]



2
different
concepts of
opacity

Article information		
	Order no. 4 g syringe	Order no. 20 x 0.3 g OptiTips [®]
Multi Opacity Refills		
Shade A1 E	7082310	7082410
Shade A2 E	7082311	7082411
Shade A2 D	7082320	7082420
Shade A3 E	7082312	7082412
Shade A3.5 E	7082313	7082413
Shade A3.5 D	7082321	7082421
Shade B1 E	7082323	7082423
Shade B2 E	7082314	7082414
Shade C2 E	7082315	7082415
Shade BL E	7082316	7082416
Shade BL D	7082322	7082422
Shade INC	7082330	7082430
Multi Opacity Kits		
Harvard PremiumFill⁺ „Multi Opacity“ Assortment Kit	7082502	
11 syringes ea. 4 g, in the shades A1 E, A2 E, B2 E, C2 E, A3 E, A3.5 E, BL E INC, A2 D, A3.5 D, BL D		
Harvard PremiumFill⁺ „Multi Opacity“ Trial Kit S	7082501	
3 syringes ea. 4 g, in the shades INC, A2 E, A2 D		
	Order no. 4 g syringe	Order no. 20 x 0.3 g OptiTips [®]
Universal Opacity Refills		
Shade A1 U	7082300	7082400
Shade A2 U	7082301	7082401
Shade A3 U	7082302	7082402
Shade A3.5 U	7082303	7082403
Shade B1 U	7082304	7082404
Shade B2 U	7082305	7082405
Shade C2 U	7082306	7082406
Universal Opacity Kits		
Harvard PremiumFill⁺ „Universal Opacity“ IntroKit S	7082500	
3 syringes ea. 4 g, in the shades A2 U, A3 U, A3.5 U		
Harvard Applier OptiTips[®]	7095200	

HARVARD[®]

Marke und Qualität seit 1892

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Phone: + 49 (0) 30/99 28 978-0
Fax: + 49 (0) 30/99 28 978-19
info@harvard-dental-international.de

11013_EN_07/18

HARVARD[®]

Harvard PremiumFill ⁺

Nano-Optimized Hybrid Composite

For **posterior** and anterior restorations.



Please rate this product.*

1. Light Cure Composites, Bondings and Liner



Harvard UltraFill

Harvard PremiumFill

Harvard UltraFill Flow

Harvard PremiumFlow

Harvard InterLock® ONE Universal Adhesive

Harvard InterLock® Mono

Harvard IonoLine

Harvard Restore

Harvard RestoreFlow

Harvard Self-Bond

Harvard Bond TE Mono

Harvard Etch

Harvard ZirconCore

Harvard Bond SE Dual

Harvard LuteCem SE

Harvard InterLock® ONE Self Cure Activator

Harvard MTA XR Flow EWT OptiCaps® & EndoDirect syringe

Harvard MTA XR Fast OptiCaps®

Harvard MTA XR OptiCaps®

2. Dual Cure Core Build-Up and Post Cementing



3. Endodontics, MTA




Harvard TEMP C&B

Harvard TEMP Glaze LC

Harvard TEMP Cem Esthetic

Harvard TEMP Cem EasyDose®


4. Temporary Materials




5.1 Zinc Oxide Luting Cements



5.2 Glass Ionomer Cements for Restorations



5.3 Glass Ionomer Cements for Luting



6. Bite Registration and Impression Materials



Harvard PremiumSIL Light/Fast, Medium/Medium Fast, Mono/Mono Fast, Heavy/Heavy Fast



Harvard TraySive



Please rate this product.*

*Great that you have tested our product. Did it meet your satisfaction? Your opinion is very important to us; please use the QR code or our website www.harvard-dental.de/feedback to provide us with your much valued feedback. Product tests are worth it. As a gesture of thanks we will send you another product sample from our range. Thank you for your time.