






HARVARD®

Harvard

POWERmer Cem



EN

Harvard POWERmer Cem

High strength, self adhesive, dual cure resin modified glass ionomer cement for luting of crowns and bridges, inlays, onlays and veneers

- >> Mild
- >> Dual cure
- >> No shrinkage
- >> **Good adhesion** to tooth, zirconium dioxide and metal
- >> **No adhesive necessary!**

- >> **Stress-bearing**
- >> **Strong like a composite cement**
- >> **High fluoride release**
- >> Insoluble
- >> Neutralises acids

- >> **Translucent**
- >> Colour stable
- >> Radiopaque

- >> HandMix or OptiCaps®

Order Details

15 g powder / 7 ml liquid, dosage spoon, mixing pad	
Universal	7071322

10 OptiCaps® ea. 0.4 g, ea. packed in aluminum pouch	
Universal	7071321

50 OptiCaps® ea. 0.4 g, ea. packed in aluminum pouch	
Universal	7071320

Harvard Applier OptiCaps®	7092000
---------------------------	---------




Setting times

	Harvard Cement NH (Harvard)	POWERmer Cem (Harvard)	RelyX Unicem Aplicap (3M)	FuJi Plus Capsule (GC)
Working time from start of mixing at 23 °C (73 °F)	2:00 min	2:30 min	2:10 min	1:40 min
Net setting time	6:00 min	2:40 min	3:20 min	2:10 min

Mechanical properties

	Harvard Cement (Harvard)	POWERmer Cem (Harvard)	RelyX Unicem Aplicap (3M)	FuJi Plus Capsule (GC)
Compressive strength: (MPa)	110	250	250	120
Flexural strength: (MPa)	10	80	80	50
E-Modul: (MPa)	15.000	10.000	9.000	7.000
Opacity: (%)	100	80	80	90
Radiopacity: (% Al)	>500	240	235	250

Fluoride release





HARVARD®

Much more
than you
expect.

Harvard Distribution Partner.

11085_EN_03/25



HARVARD®

Marke und Qualität seit 1892

Harvard Dental International GmbH
Margaretenstr. 2 - 4, 15366 Hoppegarten, Germany
Phone: + 49 (0) 30/99 28 978-0
Fax: + 49 (0) 30/99 28 978-19
info@harvard-dental-international.de



www.harvard-dental.de

Liability is excluded for all printing errors and omissions. Before using our Harvard products, the respective directions for use should be noticed in every case. All measurements are internal measurements of Harvard Dental International.