

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 25.04.2017  
 Printing date: 16.06.2017

Version: 1

### 1. Identification of the substance/mixture and of the company/undertaking

#### Product Identifier

**Trade name:** Harvard InterLock® ONE

**Use of the preparation:** Dental adhesive

**Details of the supplier of the safety data sheet:**

#### Manufacturer/supplier:

Harvard Dental International GmbH,

Margaretenstraße 2-4,

15366 Hoppegarten, Germany

info@harvard-dental-international.de

#### Emergency telephone number:

+ 49 (0) 30 99 28 978-0 (Quality Management)

### 2. Hazards identification

#### Classification of the substance or mixture:

**Classification according to Regulation (EC) No 1272/2008:**

Skin Irrit. 2, H315 Causes skin irritation

Eye Irrit. 2, H319 Causes serious eye irritation

Skin Sens1, H317 May cause allergic skin reaction

#### Information concerning particular hazards for human and environment:

The product has to be labeled due to calculation procedure for the "General Classification guideline for preparations of the EU" in the latest version.

#### Label elements

**Labeling according to Regulation (EC) 1272/2008:**

#### Hazard pictograms



#### Signal word

Warning

#### Hazard statements

H315 Causes skin irritation

H319 Causes serious eye irritation

H317 May cause allergic skin reaction

#### Precautionary statements

P280 Wear protective gloves/protective clothing/eye protection/face protection

P305+P351+P338 IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P302+P352 IF ON SKIN: Wash with plenty of soap and water.

#### Labeling according to EU guidelines:

Medical devices as defined in Directive 93/42/EEC and which are invasive or used in direct physical contact with the human body; are exempted from the provisions of Regulation (EC) No. 1272/2008 (CLP/GHS) usually if they are in the finished state and intended for the final user.

#### Other hazards

#### Results of PBT and vPvB assessment

**PBT:** not applicable

**vPvB:** not applicable

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 25.04.2017  
 Printing date: 16.06.2017

Version: 1

<b>3. <u>Composition/Information on Ingredients</u></b>						
<p><b>Chemical characterization:</b>  <b>Description:</b>          Modified methacrylates, catalysts, stabilizers in ethanol  <b>Hazardous components:</b></p> <table border="0"> <tr> <td style="padding-right: 20px;">CAS 868-77-9</td> <td style="padding-right: 20px;">2-Hydroxyethyl-methacrylate</td> <td>GHS07, H315, H319, H317</td> </tr> <tr> <td>CAS 64-17-5</td> <td>Ethanol</td> <td>GHS02, H225</td> </tr> </table>	CAS 868-77-9	2-Hydroxyethyl-methacrylate	GHS07, H315, H319, H317	CAS 64-17-5	Ethanol	GHS02, H225
CAS 868-77-9	2-Hydroxyethyl-methacrylate	GHS07, H315, H319, H317				
CAS 64-17-5	Ethanol	GHS02, H225				
<b>4. <u>First-aid measures</u></b>						
<p><b>General information:</b> Seek medical advice, if symptoms occur, which can be caused by the product  <b>after skin contact:</b> wash with plenty of water and soap  <b>after inhalation:</b> fresh air  <b>after eye contact:</b> Rinse opened eye for several minutes under running water and consult a doctor  <b>after swallowing:</b> Consult a physician.</p>						
<b>5. <u>Fire fighting measures</u></b>						
<p><b>Extinguishing media:</b>          CO<sub>2</sub>, powder, water spray.          Use extinguishing method and media suitable to surrounding fire.  <b>Special hazards arising from the substance or mixture</b> No further relevant information available  <b>Protective equipment</b> Wear self-contained respiratory protective equipment</p>						
<b>6. <u>Accidental release measures</u></b>						
<p><b>Personal precautions, protective equipment and emergency procedures:</b>          Do not inhale vapours/aerosols. Avoid substance contact. Ensure supply of fresh air in enclosed rooms  <b>Environmental precautions:</b>          Do not allow to enter drains/surface water/ ground-water.  <b>Methods for cleaning up:</b>          Absorb with liquid binding material. Clean up affected area.  <b>Reference to other sections:</b>          See Section 7 for information on safe handling          See Section 8 for information on personal protection equipment          See Section 13 for disposal information</p>						
<b>7. <u>Handling and storage</u></b>						
<p><b>Handling:</b>          This product should only be supplied or handled to dentists or dental technicians or upon their instructions.  <b>Precautions for safe handling:</b> Observe normal care for working with chemicals  <b>Information about fire and explosion protection:</b> Keep away from heat and fire. Take measure to prevent electrostatic charge  <b>Storage:</b> Store at a dry place. Store tightly closed. Storage temperature &lt;25 °C.  <b>Requirements at storerooms and containers:</b> Store only in the original containers  <b>Further information on storage conditions:</b>          Store away from food stuffs.  <b>Specific end uses:</b> No further relevant information available</p>						

# Safety Data Sheet

## according to 1907/2006/EC, Article 31

Date of issue: 25.04.2017  
 Printing date: 16.06.2017

Version: 1

### 8. Exposure controls/personal protection

**Components with limits of values to be supervised at the workplace:**

CAS-No.	name	art	value	unit
64-17-5	Ethanol	TRGS 900	500	ml/m <sup>3</sup> (ppm)
			or 960	mg/m <sup>3</sup>

**Additional information:** The lists valid during the making were used as basis

**Personal protective equipment:**

**General protective and hygienic measures:**

The usual precautionary measures are to be adhered to when handling chemicals.

Keep away from foodstuffs and beverages. Avoid contact with eyes or skin

**Hand protection**

Protective gloves.

Due to missing tests no recommendation to the glove material can be given for the product.

**Eye protection:** safety glasses

**Respiratory protection:** Required, when aerosols are generated

### 9. Physical and chemical properties

**Information on basic physical and chemical properties**

**Appearance**

<b>Form:</b>	liquid						
<b>Color:</b>	colorless						
<b>Odor:</b>	characteristic						
	<table border="1"> <thead> <tr> <th>Value</th> <th>Unit</th> <th>Method</th> </tr> </thead> <tbody> <tr> <td> </td> <td> </td> <td> </td> </tr> </tbody> </table>	Value	Unit	Method			
Value	Unit	Method					

**Safety relevant basic data**

<b>Melting point/melting range:</b>	n.a.
<b>Boiling point/boiling range:</b>	n.a.
<b>Flash point:</b>	n.a.
<b>Ignition temperature:</b>	No information available
<b>Danger of explosion:</b>	No information available
<b>Density:</b>	1.07 g/cm <sup>3</sup>
<b>Vapor pressure:</b>	n.d.
<b>pH:</b>	2 - 3
<b>Solubility in/miscibility with:</b>	
<b>Water:</b>	soluble
<b>Content of solvents:</b>	
<b>Organic solvents:</b>	ethanol

### 10. Stability and reactivity

**Thermal decomposition / conditions to be avoided**

**Possibility of hazardous reactions:** no information available

**Conditions to avoid:** Heat, direct sunlight and radical starters. The product is stable by proper use. Avoid contact with amines, peroxi compounds, polymerisation initiators, acids, bases, oxidizing agents.

**Hazardous decomposition products:** No information available.

## Safety Data Sheet

### according to 1907/2006/EC, Article 31

Date of issue: 25.04.2017  
 Printing date: 16.06.2017

Version: 1

#### **11. Toxicological information**

**Acute toxicity: (ethanol)**

LD<sub>50</sub> (oral, rat)          6200 mg/kg

**Primary irritation:**

**skin:** irritant

**eye:** irritant

**Sensitization:** sensitization by skin contact possible

**Additional toxicological information:**

Allergic reactions to methacrylates after long-term exposure are possible. When used and handled according to specifications, the product does not have any harmful effects to our experience and the information provided to us.

#### **12. Ecological information**

**General information:** Avoid transfer into environment

**Classification of water endangerment:** WGK=1 (slightly polluting substance).

#### **13. Disposal considerations**

**Product:**

**Recommendation:**

Use proper landfill disposal or incineration in accordance with local, state and federal regulations.

**Uncleaned packaging:**

**Recommendation:** Disposal must be made according to official regulations

#### **14. Transport information**

**Land transport ADR/RID and GGVS/GGVE (Germany):** no subject to transport regulations

**Sea transport IMDG:** no subject to transport regulations

**Air transport ICAO-TI and IATA-DGR:** no subject to transport regulations

#### **15. Regulatory information**

**Classification according to EC-guidelines:**

The product is a medical device according to EC-directive 93/42 EEC.

**Chemical safety assessment:**

A Chemical Safety Assessment has not been carried out.

#### **16. Other information:**

**Changes compared with the previous version:**

The above information is based on our present day knowledge and relates solely to the safety requirements of the product. The data do not signify any warranty with regards to products properties. However users of the product should satisfy themselves that the information given is sufficient and correct for their specific circumstances of use.